




BOECO BLOOD COUNTING CHAMBERS



Description	Code	Pack of
Neubauer improved, double ruling without coverglasses	BOE 01	1
Neubauer improved, double ruling with 2 coverglasses	BOE 13	1
Neubauer improved, bright-line double ruling, without coverglasses	BOE 02	1
Neubauer improved, bright-line, double ruling, with 2 coverglasses	BOE 14	1
Neubauer improved, bright-line, double ruling, with V-slash with 2 coverglasses	BOE 14V	1

We can also deliver all other types like: Neubauer, Thoma, Thoma new, Bürker, Türk, Nageotte, Fuchs-Rosenthal, Malassez.
As well as complete HAEMACYTOMETER SETS



BOECO HAEMATOLOGY PIPETTES

Description	Code	Pack of
Blood diluting pipette, acc. to Thoma, red, 1 : 11, 13 cm, blue graduation	BOE 04	10
Blood diluting pipette, acc. to Thoma, white, 1 : 11, 13 cm, blue graduation	BOE 05	10
red mouthpiece and rubber tubing	BOE 19	
white mouthpiece and rubber tubing	BOE 20	

BOECO HAEMATOLOGY PIPETTES

Description	Code	Pack of
SAHLI Haemoglobin pipette w/o funnel top, mark at 0,02 ml in 20 C°.	BOE 06	25
Haematocrit tubes acc. to WINTROBE length 110 mm, 0-100 mm in 1 mm double scale, blue graduation	BOE 17	10
Blood sedimentation pipettes acc. to WESTERGREEN, macro method length 300 mm, 0-200 mm in 1 mm white graduation	LMS 63722101	100

PASTEUR PIPETTES, SODAGLASS

Description	Code	Pack of
150 mm long	BOE 5010601	4 x 250 pcs.
230 mm long	BOE 5010602	4 x 250 pcs.


TEST TUBES

Round bottom, without rim, borosilicate glass 5.1.

Ø x L mm	Wall thickness	Code	Pack of
12 x 75	0,6 mm	BOE 03.3411.1118	250
13 x 100	0,5-0,6 mm	BOE 03.3411.1019	250
16 x 160	0,6 mm	BOE 03.3411.3020	250

Round bottom, without rim, heavy wall, with printed logo "BOECO GERMANY" and matt label, borosilicate glass 5.1.

Ø x L mm	Wall thickness	Code	Pack of
12 x 75	1,0 mm	BOE 03.3421.1126	250
13 x 100	1,0 mm	BOE 03.3421.1128	250
16 x 160	1,0±0,15 mm	BOE 03.3421.2038	250




TEST TUBES WITH THREAD

without cap, borosilicate glass 5.1.

Ø x L mm	Code	Pack of
13 x 100	BOE 155030	1000
16 x 100	BOE 155060	1000
16 x 125	BOE 155061	1000
16 x 150	BOE 155062	1000
20 x 125	BOE 155000	500
20 x 150	BOE 155001	500

SCREW CAPS, PP, BLACK, FOR TEST TUBES WITH THREAD

for tubes with Ø	Code	Pack of
13 mm	BOE 155403	1000
16 mm	BOE 155406	1000
20 mm	BOE 155409	500



VIALS, FLAT BOTTOM, WITH THREAD

without cap, borosilicate glass 5.1.

Ø x L mm, ml	clear glass	amber glass	Pack of
20 x 36, 5 ml	BOE 152035	BOE 152035B	225
23 x 38, 7,5 ml	BOE 152338	BOE 152338B	180
23 x 48, 12 ml	BOE 152348	BOE 152348B	180
23 x 72, 20 ml	BOE 152372	BOE 152372B	171

SCREW CAPS, PP, GL 18, BLUE WITH POLYETHYLENE SEALING

for tubes with Ø	Code	Pack of
GL 18	BOE 150018	100



GRADUATED PIPETTES

"BOECO GERMANY", total delivery, zero at top, from 5 ml with filter bed, amber graduation, sodaglass, Cl. AS / ISO 835

Capacity	Subdivision	Code	Pack of
0,5 ml	0,01	BOE 5116312	50
1 ml	0,01	BOE 5116314	50
2 ml	0,02	BOE 5116516	50
5 ml	0,05	BOE 5116618	50
5 ml	0,1	BOE 5116718	50
10 ml	0,1	BOE 5116719	50
20 ml	0,1	BOE 5116722	25
25 ml	0,1	BOE 5116723	25
50 ml	0,2	BOE 5116825	5
100 ml	0,2	BOE 5116831	5

VOLUMETRIC PIPETTES

"BOECO GERMANY", with one mark, amber stain graduation, sodaglass, Cl. AS / ISO 648

Capacity	Tolerance +/-ml	Code	Pack of
1 ml	0,008	BOE 5111014	50
2 ml	0,010	BOE 5111016	50
3 ml	0,010	BOE 5111017	50
4 ml	0,015	BOE 5111117	25
5 ml	0,015	BOE 5111018	25
10 ml	0,020	BOE 5111019	25
15 ml	0,030	BOE 5111021	25
20 ml	0,030	BOE 5111022	10
25 ml	0,030	BOE 5111023	10
50 ml	0,050	BOE 5111025	1
100 ml	0,080	BOE 5111030	1
200 ml	0,100	BOE 5111037	1

