


## DIMROTH CONDENSERS

borosilicate glass, olive diam. 11 mm


| Jacket length | Cone/socket size | Code         | Pack of |
|---------------|------------------|--------------|---------|
| 160 mm        | NS 14/23         | LMS 61453120 | 4       |
| 160 mm        | NS 19/26         | LMS 61453121 | 4       |
| 250 mm        | NS 19/26         | LMS 61453122 | 2       |
| 250 mm        | NS 24/29         | LMS 61453123 | 2       |
| 250 mm        | NS 29/32         | LMS 61453124 | 2       |
| 400 mm        | NS 24/29         | LMS 61453125 | 1       |
| 400 mm        | NS 29/32         | LMS 61453126 | 1       |



## JACKET COILED CONDENSERS

borosilicate glass


| Jacket length | Cone/socket size | Code         | Pack of |
|---------------|------------------|--------------|---------|
| 200 mm        | NS 29/32         | LMS 61454119 | 2       |
| 250 mm        | NS 29/32         | LMS 61454120 | 2       |
| 300 mm        | NS 29/32         | LMS 61454121 | 2       |
| 400 mm        | NS 29/32         | LMS 61454122 | 2       |



## SOXHLET EXTRACTORS

cones NS 29/32, (NS 24/40 for L version) borosilicate glass

| Vol. extractor | socket size | Code          | Pack of |
|----------------|-------------|---------------|---------|
| 30 ml          | NS 29/32    | LMS 61435100  | 10      |
| 50 ml          | NS 29/32    | LMS 61435101  | 10      |
| 100 ml         | NS 45/40    | LMS 61435102  | 10      |
| 100 ml         | NS 24/40    | LMS 61435102L | 10      |
| 150 ml         | NS 45/40    | LMS 61435103  | 10      |
| 250 ml         | NS 45/40    | LMS 61435104  | 10      |
| 250 ml         | NS 24/40    | LMS 61435104L | 10      |
| 500 ml         | NS 60/46    | LMS 61435105  | 10      |



## SOXHLET EXTRACTION APPARATUS

complete with: Soxhlet extractor, flat bottom flask NS 29/32, Dimroth condenser NS 29/32, reduction adapter

| Vol. extractor | flat bottom flask | Code          | Pack of |
|----------------|-------------------|---------------|---------|
| 30 ml          | 100 ml            | LMS 61612102  | 1       |
| 50 ml          | 250 ml            | LMS 61612103  | 1       |
| 100 ml         | 250 ml            | LMS 61612104  | 1       |
| 150 ml         | 500 ml            | LMS 61612105  | 1       |
| 250 ml         | 500 ml            | LMS 61612106  | 1       |
| 500 ml         | 1000 ml           | LMS 61612107  | 1       |
| 500 ml         | 2000 ml           | LMS 616121072 | 1       |

## WEIGHING BOTTLES

with interchangeable knob lid, matt label, low form

| Ø x Height | Code         | Pack of |
|------------|--------------|---------|
| 25 x 25 mm | LMS 61135102 | 10      |
| 35 x 30 mm | LMS 61135106 | 10      |
| 45 x 30 mm | LMS 61135109 | 10      |
| 50 x 30 mm | LMS 61135111 | 10      |
| 60 x 30 mm | LMS 61135113 | 10      |
| 80 x 30 mm | LMS 61135116 | 10      |

