

TIMER

MECHANICAL TIMER

0-60 min, sec. alarm, housing colour white/blue

Code

BOE 124



DIGITAL STOPWATCH

waterproof,
Stopwatch precision 1/100 seconds first 30 min, 1 sec. up to max 24 h,
2 splits memory
Chronograph up to 23 hours, 59 minutes, 59,99 seconds
Calendar shows day, week, month and year
Daily alarm, hourly chime and snooze function
Time shows second, minute, hour, AM/PM & 12/24 switch
With batteries (button cells, 2pcs)“

Code

BOE 260



DIGITAL TIMER

Count down 99 min., 59 sec., alarm when countdown reaches
“zero”, magnet, silver, with battery (1 pc. button cell SR 44 1,5 V.)

Code

BOE 200



DIGITAL TIMER

Count down/up timer, 99 min. 59 sec. alarm when countdown reaches
“zero”, magnet, blue, with battery (AAA, 1 pc.)

Code

BOE 210



2-LINE DIGITAL TIMER, JUMBO DISPLAY

2 independent timers in one model. Count down/Count up timer and
clock in one. Count down/up 99 hrs. 59 min., 59 sec., alarm when count-
down reaches “zero”, clip for hanging and standing, magnet,
with battery (AAA 1,5 V 1 pc.)

Code

BOE 220



3-LINE DIGITAL TIMER, JUMBO DISPLAY

Displays current time, alarm time and month & date or 3 channels timer
timing. Count down/up 99 hrs. 59 min., 59 sec., alarm when countdown
reaches “zero”, distinct alarms for each channel. Clip for hanging and stand-
ing, magnet, with battery (AAA 1,5 V 1 pc.)

Code

BOE 225

