


BOECO GLASSWARE

BEAKERS


low form, with lip and graduation, "BOECO Germany", boro 3.3.



Capacity	Code	Pack of
25 ml	BOE 5010614	10
50 ml	BOE 5010617	10
100 ml	BOE 5010624	10
150 ml	BOE 5010629	10
250 ml	BOE 5010636	10
400 ml	BOE 5010641	10
600 ml	BOE 5010648	10
800 ml	BOE 5010653	10
1000 ml	BOE 5010654	10
2000 ml	BOE 5010663	4
3000 ml	BOE 5010668	4
5000 ml	BOE 5010673	2

BEAKERS


tall form, with lip and graduation, borosilicate glass 3.3.



Capacity	Code	Pack of
50 ml	BOE 5011617	10
100 ml	BOE 5011624	10
150 ml	BOE 5011629	10
250 ml	BOE 5011636	10
400 ml	BOE 5011641	10
600 ml	BOE 5011648	10
800 ml	BOE 5011653	10
1000 ml	BOE 5011654	10
2000 ml	BOE 5011663	6
3000 ml	BOE 5011668	6

ERLENMEYER FLASKS

narrow neck flange-edged, graduated, "BOECO Germany" boro 3.3.



Capacity	Code	Pack of
25 ml	BOE 5021614	10
50 ml	BOE 5021617	10
100 ml	BOE 5021624	10
200 ml	BOE 5021632	10
250 ml	BOE 5021636	10
300 ml	BOE 5021639	10
500 ml	BOE 5021644	10
1000 ml	BOE 5021654	10
2000 ml	BOE 5021663	6
3000 ml	BOE 5021668	6
5000 ml	BOE 5021673	2

ERLENMEYER FLASKS

wide neck, flange-edged, graduated, borosilicate glass 3.3.

Capacity	Code	Pack of
25 ml	BOE 5022614	10
50 ml	BOE 5022617	10
100 ml	BOE 5022624	10
200 ml	BOE 5022632	10
250 ml	BOE 5022636	10
300 ml	BOE 5022639	10
500 ml	BOE 5022644	10
1000 ml	BOE 5022654	10
2000 ml	BOE 5022663	6



ROUND BOTTOM FLASKS

(Boiling Flasks), "BOECO Germany", boro 3.3., ISO 1773

Capacity	Code	Pack of
250 ml	BOE 5071136	10
500 ml	BOE 5071144	10



FLAT BOTTOM FLASKS

narrow neck, with beaded rim, "BOECO Germany", boro 3.3

Capacity	Code	Pack of
250 ml	BOE 5072136	10
500 ml	BOE 5072144	10



LABORATORY BOTTLES

"BOECO GERMANY" boro 3.3., graduated. Complete with blue screw cap and pouring ring of PP. Resistant to steam sterilisation up to 140 °C

Capacity	Code	Pack of
100 ml	BOE 5080245	10
250 ml	BOE 5080365	10
500 ml	BOE 5080445	10
1000 ml	BOE 5080545	10
2000 ml	BOE 5080635	10
5000 ml	BOE 5080735	6
10000 ml	BOE 5080865	1
20000 ml	BOE 5080915	1



LABORATORY BOTTLES

"BOECO GERMANY" boro 3.3., graduated. Complete with red screw cap (PBTB), PTFE protected seal and pouring ring of ETFE. Resistant to steam sterilisation up to 200 °C

Capacity	Code	Pack of
100 ml	BOE 5080241	10
250 ml	BOE 5080361	10
500 ml	BOE 5080441	10
1000 ml	BOE 5080541	10
2000 ml	BOE 5080631	10
5000 ml	BOE 5080731	6
10000 ml	BOE 5080861	1
20000 ml	BOE 5080911	1



LABORATORY BOTTLES


"BOECO GERMANY" amber, boro 3.3., graduated. Complete with blue screw cap and pouring ring of PP.
Resistant to steam sterilisation up to 140 °C



Capacity	Code	Pack of
100 ml	BOE 5086245	10
250 ml	BOE 5086365	10
500 ml	BOE 5086445	10
1000 ml	BOE 5086545	10
2000 ml	BOE 5086635	10
5000 ml	BOE 5086735	6
10000 ml	BOE 5086865	1
20000 ml	BOE 5086915	1

LABORATORY BOTTLES

"BOECO GERMANY", boro 3.3., graduated.
without screw cap and pouring ring



Capacity	Code	Pack of
100 ml	BOE 5080124	10
250 ml	BOE 5080136	10
500 ml	BOE 5080144	10
1000 ml	BOE 5080154	10
2000 ml	BOE 5080163	10
5000 ml	BOE 5080173	6
10000 ml	BOE 5080186	1
20000 ml	BOE 5080191	1




BLUE SCREW CAP (PP)

heat resistant up to 140°C
GL 45 BOE 5923928

BLUE POURING RING (PP)

heat resistant up to 140°C
GL 45 BOE 5924228



RED SCREW CAP (PBTB)


heat resistant up to 200°C
GL 45 BOE 5924028

PTFE PROTECTED SILICON SEAL

heat resistant up to 200°C
GL 45 BOE 9180001195

RED POURING RING (ETFE)

heat resistant up to 200°C
GL 45 BOE 5924428



FILTERING FLASKS


Erlenmeyer shape, SIMAX borosilicate glass
with glass hose connection, for vacuum use

Capacity	Code	Pack of
100 ml	BOE 5120124	10
250 ml	BOE 5120136	10
500 ml	BOE 5120144	8
1000 ml	BOE 5120154	8
2000 ml	BOE 5120163	6

DESICCATORS

with plastic knob and stopcock borosilicate glass 3.3.,
without porcelain plate


Capacity	Code	Pack of
100 mm	GLA 415262112	1
150 mm	GLA 415262162	1
200 mm	GLA 415262212	1
250 mm	GLA 415262262	1
300 mm	GLA 415262303	1



DESICCATORS

with plastic knob borosilicate glass 3.3., without porcelain plate


Capacity	Code	Pack of
100 mm	GLA 415262111	1
150 mm	GLA 415262161	1
200 mm	GLA 415262211	1
250 mm	GLA 415262261	1
300 mm	GLA 415262302	1



DESICCATORS

with outlet tube and stopcock NS 24/29
borosilicate glass 3.3., without porcelain plate


Capacity	Code	Pack of
100 mm	GLA 415262120	1
150 mm	GLA 415262170	1
200 mm	GLA 415262220	1
250 mm	GLA 415262270	1
300 mm	GLA 415262330	1



DESICCATOR PLATES

porcelain, glazed

Size	Code	Pack of
90 mm	GLA 9180000358	1
140 mm	GLA 9180000354	1
190 mm	GLA 9180000355	1
240 mm	GLA 9180000356	1
280 mm	GLA 9180000357	1



FUNNELS

60°, short stem, soda glass

Size	Code	Pack of
30 mm	LMS 11241020	10
35 mm	LMS 11241030	10
40 mm	LMS 11241040	10
45 mm	LMS 11241050	10
50 mm	LMS 11241060	10
55 mm	LMS 11241070	10
60 mm	LMS 11241080	10
70 mm	LMS 11241110	10
80 mm	LMS 11241160	10
100 mm	LMS 11241190	10
120 mm	LMS 11241230	10
150 mm	LMS 11241310	10
200 mm	LMS 11241380	5





REAGENT BOTTLES

narrow mouth, round shoulder, with ground-in NS Stopper

Capacity	clearglass	brownglass	Pack of
50 ml	GLA 414102050	GLA 414202050	48
100 ml	GLA 414102100	GLA 414202100	70
250 ml	GLA 414102250	GLA 414202250	30
500 ml	GLA 414102500	GLA 414202500	15
1000 ml	GLA 414102940	GLA 414202940	12
2000 ml	GLA 414102950	GLA 414202950	6



REAGENT BOTTLES

wide mouth, round shoulder, with ground-in NS Stopper

Capacity	clearglass	brownglass	Pack of
50 ml	GLA 414104050	GLA 414204050	40
100 ml	GLA 414104100	GLA 414204100	48
250 ml	GLA 414104250	GLA 414204250	30
500 ml	GLA 414104500	GLA 414204500	15
1000 ml	GLA 414104940	GLA 414204940	12
2000 ml	GLA 414104950	GLA 414204950	6



DROPPING BOTTLES

wide mouth, round shoulder, with ground-in NS Stopper

Capacity	clearglass	brownglass	Pack of
50 ml	GLA 414109050	GLA 414209050	96
100 ml	GLA 414109100	GLA 414209100	70

PETRI DISHES PETRIQ

Diameter	Code	Pcs./carton
60 x 15 mm	PET 4006015	144
70 x 15 mm	PET 4007015	144
80 x 15 mm	PET 4008015	144
90 x 15 mm	PET 4009015	72
100 x 20 mm	PET 4010010	72
100 x 15 mm	PET 4010015	72
100 x 20 mm	PET 4010020	72
120 x 20 mm	PET 4012020	72
150 x 25 mm	PET 4015025	72
180 x 30 mm	PET 4018030	6
200 x 30 mm	PET 4020030	6
200 x 50 mm	PET 4020050	6



WATCH GLASSES

fused edges, acc. to DIN 12341


Diameter	Code	Pcs./carton
50 mm	PET 3106050	100
60 mm	PET 3008060	100
70 mm	PET 3010070	100
80 mm	PET 3012080	100
90 mm	PET 3011090	80
100 mm	PET 3013100	80
120 mm	PET 3015120	60
125 mm	PET 3016125	60
150 mm	PET 3019150	60
200 mm	PET 3027200	30



STAINING DISHES

acc. to Hellendahl, with wider upper section
with cover, for 16 slides 76 x 26 mm, sodaglass


<u>I x w x h</u>	<u>Code</u>	<u>Pack of</u>
60 x 60 x 100 mm	HAV 23253111	10



STAINING DISHES

acc. to Schieferdecker,
with cover, for 10 slides 76 x 26 mm, sodaglass

<u>I x w x h</u>	<u>Code</u>	<u>Pack of</u>
90 x 70 x 40 mm	HAV 23253201	10



STAINING JAR

acc. to Coplin, with cover, for 10 slides 76 x 26 mm, sodaglass


<u>I x w x h</u>	<u>Code</u>	<u>Pack of</u>
40 / 100 mm	HAV 23253401	10



STAINING GLASS BOX

for staining tray, with cover, sodaglass


<u>I x b x h</u>	<u>Code</u>	<u>Pack of</u>
108 x 90 x 70 mm	HAV 23253423	10



STAINING TRAY

for 10 slides 76 x 26 mm, sodaglass

<u>I x w x h</u>	<u>Code</u>	<u>Pack of</u>
88 x 40 x 70 mm	HAV 23253422	10





STAINLESS STEEL HANDLE

for staining trays

<u>Code</u>	<u>Pack of</u>
HAV 23253427	10




BOECO BLOOD COUNTING CHAMBERS



Description	Code	Pack of
Neubauer improved, double ruling without coverglasses	BOE 01	1
Neubauer improved, double ruling with 2 coverglasses	BOE 13	1
Neubauer improved, bright-line double ruling, without coverglasses	BOE 02	1
Neubauer improved, bright-line, double ruling, with 2 coverglasses	BOE 14	1


We can also deliver all other types like: Neubauer, Thoma, Thoma new, Bürker, Türk, Nageotte, Fuchs-Rosenthal, Malassez.
As well as complete HAEMACYTOMETER SETS

BOECO HAEMATOLOGY PIPETTES



Description	Code	Pack of
Blood diluting Pipette, acc. to Thoma, red, 1 : 11, 13 cm, blue graduation	BOE 04	10
Blood diluting Pipette, acc. to Thoma, white, 1 : 11, 13 cm, blue graduation	BOE 05	10
red mouthpiece and rubber tubing	BOE 19	
white mouthpiece and rubber tubing	BOE 20	

BOECO HAEMATOLOGY PIPETTES



Description	Code	Pack of
SAHLI Haemoglobin pipette w/o funnel top, mark at 0,02 ml in 20 C°.	BOE 06	25
Haematocrit tubes acc. to WINTROBE length 110mm, 0-100mm in 1mm double scale, blue graduation	BOE 17	10
Blood sedimentation pipettes acc. to WESTERGREEN, macro method length 300mm, 0-200mm in 1mm white graduation	LMS 63722101	100

PASTEUR PIPETTES, SODAGLASS,



Description	Code	Pack of
150mm long	BOE 5010601	4 x 250 pcs.
230mm long	BOE 5010602	4 x 250 pcs.