



B&C HAND REFRACTOMETER

B&C HANDREFRACTOMETER


Code **Description**
B&C 30103 Saccharic Concentration 0-32 %
Resolution: 0,2 %




Code **Description**
B&C 30104 Saccharic Concentration 28-62 %
Resolution: 0,2 %





Code **Description**
B&C 30106 Saccharic Concentration 58-90 %
58-90% Brix, 38°-43°Be', 17-27% water
Resolution: 0,5 % Brix, 0,5° Be', 1% water
(Sugar content, Baume and Water content)




Code **Description**
B&C 30107 Saccharic Concentration 0-90 %
(3 ranges) 0-42%, 42-71%, 71-90%
Resolution: 0,2 %



Code **Description**
B&C 30202 Salinity 0-28 %
Resolution: 0,2 %



Code **Description**
B&C 30301 Protein in Serum 0-12g/dl. / 0,2%
Urine Gravity 1.000-1.040 / 0,002 sg



Code **Description**
B&C 32145 Portable Digital Refractometer
Range: 0-45 % Brix
Resolution: 0,1% Accuracy: ± 0,2 %
Automatic temperature compensation
Standard temperature 20°C
Sample Volume: 0,1 - 0,3 ml
Reaction time: about 3 seconds

