

PIPETTE AIDS

BOECO PIPETTE FILLER BALLS

natural rubber, red, with 3 valves

Cat. No.	Description
BOE 100	Model STANDARD, short cone for pipettes up to 25 ml
BOE 101	Model UNIVERSAL, long cone for all pipettes



BOECO PIPETTE FILLER BALLS

BOECO PIPETTE PUMPS, PP

can be operated with one hand. Glass and plastic pipettes fit easily. Precise filling or emptying is done by rotating their actuator-wheels forward or backward. For rapid emptying simply depress plunger.

Cat. No.	Description
BOE 120	0 - 2 ml, blue
BOE 121	0 - 10 ml, green
BOE 122	0 - 25 ml, red



BOECO PIPETTE PUMPS, PP

BOECO MACRO PIPETTE CONTROLLER

offers safer, more convenient filling and dispensing, with glass and plastic pipettes. With one sensitive lever filling and delivery are easily controlled. Pipetting with ease and comfort.

- ▶ Easy and tireless handling for left and right-handed people
- ▶ An especially designed intake cone allows the use of pipettes from 0.1 ml to 100 ml
- ▶ The perfectly harmonized valve unit ensures effortless liquid intake
- ▶ Precise and high sensitive set of the meniscus
- ▶ A replaceable hydrophobic membrane filter protects the interior from penetrating liquids
- ▶ Completely autoclavable

Cat. No.	Description
BOE 160	BOECO Macro Pipette controller, complete
BOE 150-1	Membrane filter, 3 µm, sterile, hydrophobic, PTFE, pack of 10
BOE 160-2	Suction bellows
BOE 160-3	Silicone adapter
BOE 160-4	Valve unit



BOECO MACRO PIPETTE CONTROLLER

BOECO ELECTRONIC MACRO PIPET. AID

is a cordless rechargeable pipetting aid which can be used with all plastic or glass pipettes between 1-100 ml, pasteur pipettes and optional also with 5 ml and 10 ml Pipettor tips. The linear speed control allows the aspiration and dispensing to be adjusted easily and precisely to suit either large or small volume pipettes. The speed can also be fine-tuned by applying varying finger pressure to the operation buttons.

- ▶ The attached stand on the unit supports unit/pipettes when not in use.
- ▶ The unit is equipped with a low battery indicator

Cat. No.	Description
BOE 9500000	BOECO Electronic Macro Pipetting Aid, complete with AC-adapter 220V, and 0,45µm hydrophobic filters
BOE 9500001	Idem with AC-adapter 110V, U.S. Plug

Accessories

BOE 712912	Replacement filter, nonsterile (5 pcs.)
BOE 711015	Nose cone and silicon adapter, 5 ml
BOE 711016	Nose cone and silicon adapter, 10 ml
BOE 985000	5 ml Pipette tip (carton of 1000 pcs.)
BOE 980320	10 ml Pipette tip (100 pcs. box)



BOECO ELECTRONIC MACRO PIPETTING AID