



BOECO CENTRIFUGE S-6

BOECO CENTRIFUGE S-6

BOECO centrifuge S-6 is widely used in physical and chemistry analysis, biochemistry, cellular and molecular biology, clinical Labs and blood donation centers

BOECO centrifuge S-6 with fixed Angle Rotor for 6 x 15 ml centrifuge tubes (17 x 100 /120 mm), 3 adapters sets, black adapters Ø 13 x 78 mm for 10 ml / 7 ml / 5 ml blood collection tubes, white adapters Ø 13 x 63 mm for 4 ml / 3 ml / 2 ml blood collection tubes, green adapters Ø 11 x 78 mm for 2,0 / 1,5 ml microtubes, Power Supply AC 100-240 V, 50-60 Hz, Manual

- » Maintenance-free brushless motor
- » LCD Display of time and speed
- » Adjustable speed from 300 rpm to 5000 rpm
- » Time, from 15 s to 99 min
- » Lid automatically opens at the end of the run to prevent sample warming and to reduce handling time
- » Two program button P1/P2 for easy access to routine programs
- » Button rpm/rcf for easy switching between RPM and RCF display
- » Automatic rotor imbalance detection

Technical details

Max RCF:	2500xg
Power supply:	External, AC 100 - 240 V / 50/60 Hz
Rotor Capacity:	6 x 15 ml (with 3 adapter sets 10/7/5ml - 4/3/2ml -1,5/2,0ml)
Speed:	From 300 rpm to 5000 rpm
Time:	15 s – 99 min. / continuous operation
Motor:	Brushless motor
Housing:	ABS, Lid: PC
Dimension:	212 x 270 x 158 mm
Weight:	4,5 kg

Code	Description
BOE 8076000	BOECO centrifuge S-6, 100-240 V, 50/60 Hz, 5000 rpm with angle rotor for 6 x 15 ml tubes and 3 adapter sets