

BOECO MICROPLATE SHAKER PSU-2T

Specially designed for medical diagnostic. Orbital shaking platform IPP-2 with a holding mechanism for 2 standard microtiter plates.

The use of direct drive and brushless motor allows continuous mixing up to 7 days and ensures reliable, trouble-free operation.

The variable speed control provides an optimally adapted speed range from 150 to 1200 rpm with an orbit of 2 mm.

Specification:

Orbit size:	2 mm
Speed range:	150 - 1200 1/min
Speed display:	scale (min/max)
Timer:	digital, 1 min - 23h 59min / non-stop
Overall size:	220 x 205 x 90 mm
Weight:	2 kg
Power:	12 V. DC 2,0A (External power adapter)

Code	Description
BOE 8052000	Microtiter Plate Shaker PSU-2T, incl. platform IPP-2 AC Adapter 100-240 V with Multi plug (EU, US, UK, AU)
BOE 8052002	Platform IPP-4 for 4 Microtiter Plates



PSU-2T



PLATFORM IPP-4

BOECO MICROPLATE INCUBATOR SHAKER PST-60 HL

The PST-60 HL shakers provide orbital shaking of the samples and have a heating platform for 2 x 96 or 4 x 96 wells microplates and an upper heating lid. The display indicates simultaneous set and actual time, temperature and shaking speed. The PST-60 HL are ideal for a variety of cytochemistry, immunochemistry, biochemistry and research applications.

Specification:	PST-60 HL	PST-60 HL 4
Temperature regulation range:	25°C 60°C	
Temperature setting resolution:	0,1°C	
Temp. Stability:	± 0,1 °C	
Temp. Uniformity over the platform:	± 0,25 °C	
Time of thermoblock heating from RT till 37°C:	12 min.	
Orbit :	2 mm	
Max. height of micotest plate:	18 mm	
Speed regulation:	250 –1200 rpm / increment 10 rpm	
Timer with sound signal:	1 min - 96 hrs, increment 1 min	
Power (external) :	DC 12 V. 5 A	
Platform dimensions:	250 x 150 mm	210 x 290 mm
Number of microplates:	2	4
Size :	270x260x125 mm	380x390x140 mm
Weight incl. power supply :	6,1 kg	9,1 kg

Code	Description
BOE 8072000	Microplate Incubator Shaker PST-60 HL, for 2 plates, AC 100-240 V, 50/60 Hz
BOE 8072040	Microplate Incubator Shaker PST-60 HL 4, for 4 plates, AC 100V-240 V, 50/60 Hz



PST-60 HL



PST-60 HL 4

BOECO

The family-owned company Boeckel + Co. (GmbH + Co.) KG was founded in 1929 and now exists in the fourth generation.

To our customers on all 5 continents we are pleased to introduce our fifth general BOECO catalog with our extended and updated complete line of laboratory equipment and consumables.

We offer high quality at competitive prices by a convenient and cost-effective one-stop supply.

CONTACT

Boeckel + Co. (GmbH + Co.) KG
Roedingsmarkt 33
20459 Hamburg
Germany

Phone: + 49 40 325627-0
Fax: + 49 40 325627-41
E-mail: info@boeco.com
Internet: www.boeco.com

V.A.T. ID-No: DE 118903538

BANK

EURO €

Deutsche Bank AG, Hamburg,
SWIFT: DEUTDEHHXXX
IBAN- No. DE 65 2007 0000 0010 5601 00

USD \$

Deutsche Bank AG, Hamburg,
SWIFT: DEUTDEHHXXX
IBAN- No. DE 27 2007 0000 0010 5601 05