




LABORATORY GLASSWARE AND PORCELAIN

B **BOECO**
Germany

BOECO GLASSWARE

LABORATORY BOTTLES

"BOECO GERMANY" boro 3.3., graduated. Acc. to ISO 4796-1,
Complete with **blue** screw cap and pouring ring of PP
Resistant to steam sterilization up to 140 °C, Retrace code



Capacity	GL	Code	Pack of
50 ml	32	BOE 5080175	10
100 ml	45	BOE 5080245	10
250 ml	45	BOE 5080365	10
500 ml	45	BOE 5080445	10
1000 ml	45	BOE 5080545	10
2000 ml	45	BOE 5080635	10
5000 ml	45	BOE 5080735	6
10000 ml	45	BOE 5080865	1
20000 ml	45	BOE 5080915	1

LABORATORY BOTTLES

"BOECO GERMANY" boro 3.3., graduated. Acc. to ISO 4796-1,
Complete with **red** screw cap (PBTB), PTFE protected seal and pouring ring
of ETFE. Resistant to steam sterilization up to 200 °C, Retrace code



Capacity	GL	Code	Pack of
100 ml	45	BOE 5080241	10
250 ml	45	BOE 5080361	10
500 ml	45	BOE 5080441	10
1000 ml	45	BOE 5080541	10
2000 ml	45	BOE 5080631	10
5000 ml	45	BOE 5080731	6
10000 ml	45	BOE 5080861	1
20000 ml	45	BOE 5080911	1

LABORATORY BOTTLES


"BOECO GERMANY" boro 3.3., amber, graduated. Acc. to ISO 4796-1,
Complete with **blue** screw cap and pouring ring of PP.
Resistant to steam sterilization up to 140 °C, Retrace code



Capacity	GL	Code	Pack of
50 ml	32	BOE 5086175	10
100 ml	45	BOE 5086245	10
250 ml	45	BOE 5086365	10
500 ml	45	BOE 5086445	10
1000 ml	45	BOE 5086545	10
2000 ml	45	BOE 5086635	10
5000 ml	45	BOE 5086735	6
10000 ml	45	BOE 5086865	1
20000 ml	45	BOE 5086915	1

LABORATORY BOTTLES, WIDE NECK

"BOECO GERMANY" boro 3.3., graduated, square.
Complete with **blue** GL 80 screw cap and pouring ring of PP.
Resistant to steam sterilization up to 140 °C, Retrace code




Capacity	GL	Code	Pack of
500 ml	80	BOE 5098445	10
1000 ml	80	BOE 5098545	10
2000 ml	80	BOE 5098635	10
5000 ml	80	BOE 5098735	6

LABORATORY BOTTLES, WIDE NECK

"BOECO GERMANY" boro 3.3., amber, graduated, square, Retrace code
Complete with **blue** GL 80 screw cap and pouring ring of PP.
Resistant to steam sterilization up to 140 °C

Capacity	GL	Code	Pack of
500 ml	80	BOE 5096445	10
1000 ml	80	BOE 5096545	10
2000 ml	80	BOE 5096635	10
5000 ml	80	BOE 5096735	6



SCREW CAPS

GL	Color, material	Heat resistant	Code	Pack of
32	Blue, PP	up to 140°C	BOE 5923919	10
45	Blue, PP	up to 140°C	BOE 5923928	10
45	Red, PBTB	up to 200°C	BOE 5924028	10
80	Blue, PP	up to 140°C	BOE 5923916	10



POURING/OUTLET RINGS,

GL	Color, material	Heat resistant	Code	Pack of
32	Blue, PP	up to 140°C	BOE 5924219	10
45	Blue, PP	up to 140°C	BOE 5924228	10
45	Red, ETFE	up to 200°C	BOE 5924428	10
80	Blue, PP	up to 140°C	BOE 5924207	10



PTFE COATED SILICON GASKET,

Diameter	Heat resistant	Code	Pack of
43 mm	up to 200°C	BOE 9180001239	10

BEAKERS

low form, with lip and graduation, "BOECO Germany", boro 3.3., ISO 3819


Capacity	Code	Pack of
25 ml	BOE 5010614	10
50 ml	BOE 5010617	10
100 ml	BOE 5010624	10
150 ml	BOE 5010629	10
250 ml	BOE 5010636	10
400 ml	BOE 5010641	10
600 ml	BOE 5010648	10
800 ml	BOE 5010653	10
1000 ml	BOE 5010654	10
2000 ml	BOE 5010663	4
3000 ml	BOE 5010668	4
5000 ml	BOE 5010673	2



BEAKERS


tall form, with lip and graduation, SIMAX, borosilicate glass 3.3., ISO 3819

Capacity	Code	Pack of
50 ml	BOE 5011617	10
100 ml	BOE 5011624	10
150 ml	BOE 5011629	10
250 ml	BOE 5011636	10
400 ml	BOE 5011641	10
600 ml	BOE 5011648	10
800 ml	BOE 5011653	10
1000 ml	BOE 5011654	10
2000 ml	BOE 5011663	6
3000 ml	BOE 5011668	6



ERLENMEYER FLASKS


narrow neck, flange-edged, graduated, "BOECO Germany" boro 3.3.



Capacity	Code	Pack of
25 ml	BOE 5021614	10
50 ml	BOE 5021617	10
100 ml	BOE 5021624	10
250 ml	BOE 5021636	10
500 ml	BOE 5021644	10
1000 ml	BOE 5021654	10
2000 ml	BOE 5021663	6
3000 ml	BOE 5021668	6
5000 ml	BOE 5021673	2

ERLENMEYER FLASKS

wide neck, flange-edged, graduated, SIMAX borosilicate glass 3.3.
ISO 24450



Capacity	Code	Pack of
25 ml	BOE 5022614	10
50 ml	BOE 5022617	10
100 ml	BOE 5022624	10
250 ml	BOE 5022636	10
500 ml	BOE 5022644	10
1000 ml	BOE 5022654	10
2000 ml	BOE 5022663	6

ROUND BOTTOM FLASKS


(Boiling Flasks), "BOECO Germany" borosilicate glass 3.3.



Capacity	Code	Pack of
250 ml	BOE 5072136	10
500 ml	BOE 5072144	10

FLAT BOTTOM FLASKS


narrow neck, with beaded rim, SIMAX borosilicate glass 3.3.



Capacity	Code	Pack of
250 ml	BOE 5071136	10
500 ml	BOE 5071144	10

FILTERING FLASKS

Erlenmeyer shape, SIMAX borosilicate glass 3.3.
with glass hose connection, for vacuum use




Capacity	Code	Pack of
100 ml	BOE 5120124	10
250 ml	BOE 5120136	10
500 ml	BOE 5120144	8
1000 ml	BOE 5120154	8
2000 ml	BOE 5120163	6

DESICCATORS

with plastic knob and stopcock, borosilicate glass 3.3.
without porcelain plate


Nom. size	Code	Pack of
100 mm	GLA 415262112	1
150 mm	GLA 415262162	1
200 mm	GLA 415262212	1
250 mm	GLA 415262262	1
300 mm	GLA 415262303	1



DESICCATORS

with plastic knob, borosilicate glass 3.3., without porcelain plate


Nom. size	Code	Pack of
100 mm	GLA 415262111	1
150 mm	GLA 415262161	1
200 mm	GLA 415262211	1
250 mm	GLA 415262261	1
300 mm	GLA 415262302	1



DESICCATORS

with outlet tube and stopcock NS 24/29,
borosilicate glass 3.3., without porcelain plate


Nom. size	Code	Pack of
100 mm	GLA 415262120	1
150 mm	GLA 415262170	1
200 mm	GLA 415262220	1
250 mm	GLA 415262270	1
300 mm	GLA 415262330	1



DESICCATOR PLATES

porcelain, glazed

Ø	Code	Pack of
90 mm	GLA 9180000358	1
140 mm	GLA 9180000354	1
190 mm	GLA 9180003405	1
240 mm	GLA 9180003398	1
280 mm	GLA 9180003410	1



FUNNELS

60°, short stem, borosilicate glass 3.3.

Ø	Code	Pack of
30 mm	LMS 11241020	10
35 mm	LMS 11241030	10
40 mm	LMS 11241040	10
45 mm	LMS 11241050	10
50 mm	LMS 11241060	10
55 mm	LMS 11241070	10
60 mm	LMS 11241080	10
70 mm	LMS 11241110	10
80 mm	LMS 11241160	10
100 mm	LMS 11241190	10
120 mm	LMS 11241230	10
150 mm	LMS 11241310	10
200 mm	LMS 11241380	5





REAGENT BOTTLES

narrow mouth, round shoulder, with ground-in NS Stopper, sodaglass

Capacity	clear glass	brown glass	Pack of
50 ml	GLA 414102050	GLA 414202050	10
100 ml	GLA 414102100	GLA 414202100	10
250 ml	GLA 414102250	GLA 414202250	10
500 ml	GLA 414102500	GLA 414202500	10
1000 ml	GLA 414102940	GLA 414202940	10
2000 ml	GLA 414102950	GLA 414202950	6



REAGENT BOTTLES

wide mouth, round shoulder, with ground-in NS Stopper, sodaglass

Capacity	clear glass	brown glass	Pack of
50 ml	GLA 414104050	GLA 414204050	40
100 ml	GLA 414104100	GLA 414204100	10
250 ml	GLA 414104250	GLA 414204250	10
500 ml	GLA 414104500	GLA 414204500	10
1000 ml	GLA 414104940	GLA 414204940	10
2000 ml	GLA 414104950	GLA 414204950	6



DROPPING BOTTLES

wide mouth, round shoulder, with ground-in NS Stopper, sodaglass

Capacity	clear glass	brown glass	Pack of
50 ml	GLA 414109050	GLA 414209050	96
100 ml	GLA 414109100	GLA 414209100	10



BOECO PETRI DISHES

Ø x height	Code	pcs./carton
60 x 15 mm	BOE 4006015	144
70 x 15 mm	BOE 4007015	144
80 x 15 mm	BOE 4008015	144
90 x 15 mm	BOE 4009015	72
100 x 10 mm	BOE 4010010	72
100 x 15 mm	BOE 4010015	72
100 x 20 mm	BOE 4010020	72
120 x 20 mm	BOE 4012020	72
150 x 25 mm	BOE 4015025	72
180 x 30 mm	BOE 4018030	6
200 x 30 mm	BOE 4020030	6
200 x 50 mm	BOE 4020050	6



WATCH GLASSES


fused edges, acc. to DIN 12341

Ø	Code	pcs./carton
50 mm	BOE 3106050	100
60 mm	BOE 3008060	100
70 mm	BOE 3010070	100
80 mm	BOE 3012080	100
90 mm	BOE 3011090	80
100 mm	BOE 3013100	80
120 mm	BOE 3015120	60
125 mm	BOE 3016125	60
150 mm	BOE 3019150	60
200 mm	BOE 3027200	30

STAINING DISHES

acc. to Hellendahl, with wider upper section,
with cover, for 16 slides 76 x 26 mm, sodaglass


l x w x h	Code	Pack of
60 x 60 x 100 mm	HAV 23253111	10



STAINING DISHES

acc. to Schieferdecker,
with cover, for 10 slides 76 x 26 mm, sodaglass


l x w x h	Code	Pack of
90 x 70 x 40 mm	HAV 23253201	10



STAINING JAR

acc. to Coplin, with cover, for 10 slides 76 x 26 mm, sodaglass


w x h	Code	Pack of
40 / 100 mm	HAV 23253401	10



STAINING GLASS BOX

for staining tray, with cover, sodaglass


l x w x h	Code	Pack of
108 x 90 x 70 mm	HAV 23253423	10



STAINING TRAY

for 10 slides 76 x 26 mm, sodaglass

l x w x h	Code	Pack of
88 x 40 x 70 mm	HAV 23253422	10





STAINLESS STEEL HANDLE

for staining trays

Code	Pack of
HAV 23253427	10



BOECO BLOOD COUNTING CHAMBERS



Description	Code	Pack of
Neubauer improved, double ruling without coverglasses	BOE 01	1
Neubauer improved, double ruling with 2 coverglasses	BOE 13	1
Neubauer improved, bright-line double ruling, without coverglasses	BOE 02	1
Neubauer improved, bright-line, double ruling, with 2 coverglasses	BOE 14	1
Neubauer improved, bright-line, double ruling, with V-slash with 2 coverglasses	BOE 14V	1

We can also deliver all other types like: Neubauer, Thoma, Thoma new, Bürker, Türk, Nageotte, Fuchs-Rosenthal, Malassez.
As well as complete HAEMACYTOMETER SETS

PASTEUR PIPETTES, SODAGLASS

Description	Code	Pack of
150 mm long	BOE 5010601	4 x 250 pcs.
230 mm long	BOE 5010602	4 x 250 pcs.


TEST TUBES

Round bottom, without rim, borosilicate glass 5.1.

Ø x L mm	Wall thickness	Code	Pack of
12 x 75	0,6 mm	BOE 03.3411.1118	250
13 x 100	0,5-0,6 mm	BOE 03.3411.1019	250
16 x 160	0,6 mm	BOE 03.3411.3020	250

Round bottom, without rim, heavy wall, with printed logo "BOECO GERMANY" and matt label, borosilicate glass 5.1.

Ø x L mm	Wall thickness	Code	Pack of
12 x 75	1,0 mm	BOE 03.3421.1126	250
13 x 100	1,0 mm	BOE 03.3421.1128	250
16 x 160	1,0±0,15 mm	BOE 03.3421.2038	250




TEST TUBES WITH THREAD

without cap, borosilicate glass 5.1.

Ø x L mm	Code	Pack of
13 x 100	BOE 155030	1000
16 x 100	BOE 155060	1000
16 x 125	BOE 155061	1000
16 x 150	BOE 155062	1000
20 x 125	BOE 155000	500
20 x 150	BOE 155001	500

SCREW CAPS, PP, BLACK, FOR TEST TUBES WITH THREAD

for tubes with Ø	Code	Pack of
13 mm	BOE 155403	1000
16 mm	BOE 155406	1000
20 mm	BOE 155409	500



VIALS, FLAT BOTTOM, WITH THREAD

without cap, borosilicate glass 5.1.

Ø x L mm	ml	clear glass	amber glass	Pack of
20 x 36	5 ml	BOE 152035	BOE 152035B	225
23 x 38	7,5 ml	BOE 152338	BOE 152338B	180
23 x 48	12 ml	BOE 152348	BOE 152348B	180
23 x 72	20 ml	BOE 152372	BOE 152372B	171

SCREW CAPS, PP, GL 18, BLUE WITH POLYETHYLENE SEALING

for tubes with Ø	Code	Pack of
GL 18	BOE 150018	100



GRADUATED PIPETTES

“BOECO GERMANY”, total delivery, zero at top, from 5 ml with filter bed, amber graduation, sodaglass, Cl. AS / ISO 835

Capacity	Subdivision	Code	Pack of
0,5 ml	0,01	BOE 5116312	50
1 ml	0,01	BOE 5116314	50
2 ml	0,02	BOE 5116516	50
5 ml	0,05	BOE 5116618	50
5 ml	0,1	BOE 5116718	50
10 ml	0,1	BOE 5116719	50
20 ml	0,1	BOE 5116722	25
25 ml	0,1	BOE 5116723	25
50 ml	0,2	BOE 5116825	5
100 ml	0,2	BOE 5116831	5



VOLUMETRIC PIPETTES

“BOECO GERMANY”, with one mark, amber stain graduation, sodaglass, Cl. AS / ISO 648

Capacity	Tolerance +/-ml	Code	Pack of
1 ml	0,008	BOE 5111014	50
2 ml	0,010	BOE 5111016	50
3 ml	0,010	BOE 5111017	50
4 ml	0,015	BOE 5111117	25
5 ml	0,015	BOE 5111018	25
10 ml	0,020	BOE 5111019	25
15 ml	0,030	BOE 5111021	25
20 ml	0,030	BOE 5111022	10
25 ml	0,030	BOE 5111023	10
50 ml	0,050	BOE 5111025	1
100 ml	0,080	BOE 5111030	1
200 ml	0,100	BOE 5111037	1





BURETTES


"BOECO GERMANY" acc. to Mohr, with straight interchangeable glass stopcock borosilicate glass, CL. AS, ISO 385

Capacity	Subdiv. ml	Code	Pack of
10 ml	0,05	BOE 5156619	1
10 ml	0,1	BOE 5156719	1
25 ml	0,05	BOE 5156623	1
25 ml	0,1	BOE 5156723	1
50 ml	0,1	BOE 5156725	1
100 ml	0,2	BOE 5156831	1

SCHELLBACH BURETTES

"BOECO GERMANY" acc. to Mohr, with straight interchangeable glass stopcock borosilicate glass, CL. AS, ISO 385

Capacity	Subdiv. ml	Code	Pack of
10 ml	0,05	BOE 5176619	1
10 ml	0,1	BOE 5176719	1
25 ml	0,05	BOE 5176623	1
25 ml	0,1	BOE 5176723	1
50 ml	0,1	BOE 5176725	1
100 ml	0,2	BOE 5176831	1



BURETTES


"BOECO GERMANY" acc. to Mohr, with straight PTFE stopcock borosilicate glass, CL. AS, ISO 385

Capacity	Subdiv. ml	Code	Pack of
10 ml	0,05	BOE 5256619	1
10 ml	0,1	BOE 5256719	1
25 ml	0,05	BOE 5256623	1
25 ml	0,1	BOE 5256723	1
50 ml	0,1	BOE 5256725	1
100 ml	0,2	BOE 5256831	1

SCHELLBACH BURETTES

"BOECO GERMANY" acc. to Mohr, with straight PTFE stopcock borosilicate glass, CL. AS, ISO 385

Capacity	Subdiv. ml	Code	Pack of
10 ml	0,05	BOE 5276619	1
10 ml	0,1	BOE 5276719	1
25 ml	0,05	BOE 5276623	1
25 ml	0,1	BOE 5276723	1
50 ml	0,1	BOE 5276725	1
100 ml	0,2	BOE 5276831	1



BURETTES

"BOECO GERMANY" acc. to Mohr, with straight PTFE spindle stopcock borosilicate glass, CL. AS, ISO 385

Capacity	Subdiv. ml	Code	Pack of
10 ml	0,05	BOE 5458619	1
10 ml	0,1	BOE 5458719	1
25 ml	0,05	BOE 5458623	1
25 ml	0,1	BOE 5458723	1
50 ml	0,1	BOE 5458725	1
100 ml	0,2	BOE 5458831	1

SCHELLBACH BURETTES

"BOECO GERMANY" acc. to Mohr, with straight PTFE spindle stopcock stopcock borosilicate glass, CL. AS, ISO 385

Capacity	Subdiv. ml	Code	Pack of
10 ml	0,05	BOE 5478619	1
10 ml	0,1	BOE 5478719	1
25 ml	0,05	BOE 5478623	1
25 ml	0,1	BOE 5478723	1
50 ml	0,1	BOE 5478725	1
100 ml	0,2	BOE 5478831	1

VOLUMETRIC PIPETTES

with one mark, amber stain graduation,
sodaglass, DIN EN ISO 648


Capacity	Tolerance +/-ml	Cl. AS	Pack of
0,5 ml	0,010/0,005	LMS 63331421	50
1 ml	0,015/0,007	LMS 63331422	50
2 ml	0,020/0,010	LMS 63331423	30
3 ml	0,020/0,015	LMS 63331424	30
4 ml	0,030/0,015	LMS 63331425	15
5 ml	0,030/0,015	LMS 63331426	10
6 ml	0,030/0,020	LMS 63331427	10
10 ml	0,040/0,025	LMS 63331428	10
15 ml	0,060/0,030	LMS 63331429	20
20 ml	0,060/0,030	LMS 63331430	20
25 ml	0,060/0,030	LMS 63331431	10
30 ml	0,100/0,050	LMS 63331432	10
40 ml	0,100/0,050	LMS 63331433	10
50 ml	0,100/0,050	LMS 63331434	10
100 ml	0,150/0,080	LMS 63331435	6



GRADUATED PIPETTES

total delivery, zero at top, from 5 ml with filter bed, amber graduation,
sodaglass, DIN EN ISO 835

Capacity	Subdivision	Cl. B	Cl. AS	Pack of
0,1 ml	0,01	LMS 63336151	LMS 63336161	100
0,2 ml	0,02	LMS 63336152	LMS 63336162	100
0,5 ml	0,01	LMS 63336153	LMS 63336163	100
1 ml	0,01	LMS 63336154	LMS 63336164	100
2 ml	0,02	LMS 63336155	LMS 63336165	100
5 ml	0,05	LMS 63336156	LMS 63336166	100
5 ml	0,1	LMS 63336157	LMS 63336167	100
10 ml	0,1	LMS 63336158	LMS 63336168	100
20 ml	0,1	LMS 63336159	LMS 63336169	100
25 ml	0,1	LMS 63336160	LMS 63336170	100
50 ml	0,2	LMS 63336181	LMS 63336182	100



PYCROMETER ACC. TO GAY LUSSAC


Capacity	uncalibrated	calibrated to Contain (In)	Pack of
1 ml	LMS 61372100	LMS 61372110	2
2 ml	LMS 61372101	LMS 61372111	2
5 ml	LMS 61372102	LMS 61372112	2
10 ml	LMS 61372103	LMS 61372113	2
25 ml	LMS 61372104	LMS 61372114	2
50 ml	LMS 61372105	LMS 61372115	2
100 ml	LMS 61372106	LMS 61372116	2




PYCROMETER

acc. to Gay Lussac with lateral capillary tube and thermometer, calibrated to contain (In)

Capacity	Code	Pack of
10 ml	LMS 61372313	2
25 ml	LMS 61372314	2
50 ml	LMS 61372315	2
100 ml	LMS 61372316	2





MEASURING CYLINDER

tall form, hexagonal glass base, Cl. B,
borosilicate glass, acc. to DIN EN ISO 4788


Capacity	Subdiv. ml	blue graduation	amber stain graduation	Pack of
5 ml	0,1	LMS 61321201	LMS 61321231	10
10 ml	0,2	LMS 61321202	LMS 61321232	10
25 ml	0,5	LMS 61321203	LMS 61321233	10
50 ml	1	LMS 61321204	LMS 61321234	10
100 ml	1	LMS 61321205	LMS 61321235	10
250 ml	2	LMS 61321206	LMS 61321236	10
500 ml	5	LMS 61321207	LMS 61321237	10
1000 ml	10	LMS 61321208	LMS 61321238	10
2000 ml	20	LMS 61321209	LMS 61321239	10



MEASURING CYLINDER

tall form, hexagonal glass base, Cl. A,
borosilicate glass, acc. to DIN EN ISO 4788


Capacity	Subdiv. ml	blue graduation	amber stain graduation	Pack of
5 ml	0,1	LMS 61321021	LMS 61321031	10
10 ml	0,2	LMS 61321022	LMS 61321032	10
25 ml	0,5	LMS 61321023	LMS 61321033	10
50 ml	1	LMS 61321024	LMS 61321034	10
100 ml	1	LMS 61321025	LMS 61321035	10
250 ml	2	LMS 61321026	LMS 61321036	10
500 ml	5	LMS 61321027	LMS 61321037	10
1000 ml	10	LMS 61321028	LMS 61321038	10
2000 ml	20	LMS 61321029	LMS 61321039	10



MEASURING CYLINDER

tall form, hexagonal plastic base, Cl. B,
sodaglass, acc. to DIN EN ISO 4788

Capacity	Subdiv. ml	blue graduation	amber stain graduation	Pack of
10 ml	0,2	LMS 63321302	LMS 63321332	10
25 ml	0,5	LMS 63321303	LMS 63321333	10
50 ml	1	LMS 63321304	LMS 63321334	10
100 ml	1	LMS 63321305	LMS 63321335	10
250 ml	2	LMS 63321306	LMS 63321336	10
500 ml	5	LMS 63321307	LMS 63321337	10
1000 ml	10	LMS 63321308	LMS 63321338	10
2000 ml	20	LMS 63321309	LMS 63321339	10



MEASURING CYLINDER

tall form, hexagonal plastic base, Cl. A,
sodaglass, acc. to DIN EN ISO 4788

Capacity	Subdiv. ml	blue graduation	amber stain graduation	Pack of
10 ml	0,2	LMS 63321312	LMS 63321322	10
25 ml	0,5	LMS 63321313	LMS 63321323	10
50 ml	1	LMS 63321314	LMS 63321324	10
100 ml	1	LMS 63321315	LMS 63321325	10
250 ml	2	LMS 63321316	LMS 63321326	10
500 ml	5	LMS 63321317	LMS 63321327	10
1000 ml	10	LMS 63321318	LMS 63321328	10
2000 ml	20	LMS 63321319	LMS 63321329	10

MIXING CYLINDER

hexagonal glass base, plastic stopper,
blue graduation, borosilicate glass, acc. to DIN EN ISO 4788

Capacity	Neck size	Cl. B	Cl. A	Pack of
10 ml	NS 10/19	LMS 61322832	LMS 61322822	10
25 ml	NS 14/23	LMS 61322833	LMS 61322823	10
50 ml	NS 19/16	LMS 61322834	LMS 61322824	10
100 ml	NS 24/29	LMS 61322835	LMS 61322825	10
250 ml	NS 29/32	LMS 61322836	LMS 61322826	10
500 ml	NS 34/35	LMS 61322838	LMS 61322828	10
1000 ml	NS 45/40	LMS 61322839	LMS 61322829	10
2000 ml	NS 45/40	LMS 61322840	LMS 61322830	10





IMAGE WITH AMBER
GRADUATION

MIXING CYLINDER

hexagonal glass base, glass stopper,
blue graduation, borosilicate glass, acc. to DIN EN ISO 4788

Capacity	Neck size	Cl. B	Cl. A	Pack of
10 ml	NS 10/19	LMS 61322812	LMS 61322852	10
25 ml	NS 14/23	LMS 61322813	LMS 61322853	10
50 ml	NS 19/16	LMS 61322814	LMS 61322854	10
100 ml	NS 24/29	LMS 61322815	LMS 61322855	10
250 ml	NS 29/32	LMS 61322816	LMS 61322856	10
500 ml	NS 34/35	LMS 61322818	LMS 61322858	10
1000 ml	NS 45/40	LMS 61322819	LMS 61322859	10
2000 ml	NS 45/40	LMS 61322820	LMS 61322860	10




SEDIMENTATION CONES

acc. to Imhoff graduation till 100 ml, ringmarks at 1000 ml
borosilicate glass, acc. to DIN 12672

Capacity	without stopcock	with stopcock	Pack of
1000 ml	LMS 61329102	LMS 61329107	2

graduation till 1000 ml


Capacity	without stopcock	with stopcock	Pack of
1000 ml	LMS 61329114	LMS 61329118	2



COLOUR COMPARISON TUBES

acc. to Nessler, sodaglass, 2 ring marks at 50 and 100 ml

Capacity	type	Code	Pack of
100 ml	low form	LMS 63324203	10
100 ml	tall form	LMS 63324103	10





STANDARD STOPPER



OCTAGONAL STOPPER



VOLUMETRIC FLASKS

Cl. A with plastic stopper, borosilicate glass, acc. to DIN EN ISO 1042

Capacity	Neck size	standard stopper	octagonal stopper	Pack of
5 ml	NS 7/16	LMS 61319620	LMS 61319520DP	10
10 ml	NS 7/16	LMS 61319621	LMS 61319521DP	10
25 ml	NS 10/19	LMS 61319622	LMS 61319522DP	10
50 ml	NS 12/21	LMS 61319623	LMS 61319523DP	10
100 ml	NS 12/21	LMS 61319624	LMS 61319524DP	10
200 ml	NS 14/23	LMS 61319625	LMS 61319525DP	10
250 ml	NS 14/23	LMS 61319626	LMS 61319526DP	10
500 ml	NS 19/26	LMS 61319627	LMS 61319527DP	10
1000 ml	NS 24/29	LMS 61319628	LMS 61319528DP	10
2000 ml	NS 29/32	LMS 61319629	LMS 61319529DP	5
5000 ml	NS 34/35	LMS 61319630	LMS 61319530DP	1

VOLUMETRIC FLASKS

Cl. A, with glass stopper, borosilicate glass, acc. to DIN EN ISO 1042

Capacity	Neck size	glass stopper	Pack of
5 ml	NS 7/16	LMS 61319520	10
10 ml	NS 7/16	LMS 61319521	10
25 ml	NS 10/19	LMS 61319522	10
50 ml	NS 12/21	LMS 61319523	10
100 ml	NS 12/21	LMS 61319524	10
200 ml	NS 14/23	LMS 61319525	10
250 ml	NS 14/23	LMS 61319526	10
500 ml	NS 19/26	LMS 61319527	10
1000 ml	NS 24/29	LMS 61319528	10
2000 ml	NS 29/32	LMS 61319529	5
5000 ml	NS 29/32	LMS 61319530	1

VOLUMETRIC FLASKS

Cl. A, with plastic stopper, amber borosilicate glass white graduation, acc. to DIN EN ISO 1042

Capacity	Neck size	standard stopper	octagonal stopper	Pack of
5 ml	NS 7/16	LMS 61319960	LMS 61319960DP	10
10 ml	NS 7/16	LMS 61319961	LMS 61319961DP	10
25 ml	NS 10/19	LMS 61319962	LMS 61319962DP	10
50 ml	NS 12/21	LMS 61319963	LMS 61319963DP	10
100 ml	NS 12/21	LMS 61319964	LMS 61319964DP	10
200 ml	NS 14/23	LMS 61319965	LMS 61319965DP	10
250 ml	NS 14/23	LMS 61319966	LMS 61319966DP	10
500 ml	NS 19/26	LMS 61319967	LMS 61319967DP	10
1000 ml	NS 24/29	LMS 61319968	LMS 61319968DP	10
2000 ml	NS 29/32	LMS 61319969	LMS 61319969DP	5

VOLUMETRIC FLASKS

Cl. A, with glass stopper, amber borosilicate glass white graduation, acc. to DIN EN ISO 1042

Capacity	Neck size	glass stopper	Pack of
5 ml	NS 7/16	LMS 61319820	10
10 ml	NS 7/16	LMS 61319821	10
25 ml	NS 10/19	LMS 61319822	10
50 ml	NS 12/21	LMS 61319823	10
100 ml	NS 12/21	LMS 61319824	10
200 ml	NS 14/23	LMS 61319825	10
250 ml	NS 14/23	LMS 61319826	10
500 ml	NS 19/26	LMS 61319827	10
1000 ml	NS 24/29	LMS 61319828	10
2000 ml	NS 29/32	LMS 61319829	5



VOLUMETRIC FLASKS

for sugar testing, with 2 marks and rim
borosilicate glass, calibrated "to contain" and marked at 20°C

Capacity	Cl. A	Pack of
50/55 ml	LMS 61319121	10
100/110 ml	LMS 61319122	10



VOLUMETRIC FLASKS

for sugar testing acc. to Kohlrausch with rim
borosilicate glass, calibrated "to contain" and marked at 20°C


Capacity	Cl. A	Pack of
100 ml	LMS 61313125	10
200 ml	LMS 61313126	10



SEPARATING FUNNELS

conical form (pear type) acc. to DIN ISO 4800, with plastic stopper,
graduation, borosilicate glass

Capacity	Neck size	PTFE stopcock key	glass stopcock key	Pack of
50 ml	NS 19/26	LMS 61165601	LMS 61165611	10
100 ml	NS 19/26	LMS 61165602	LMS 61165612	10
250 ml	NS 29/32	LMS 61165603	LMS 61165613	10
500 ml	NS 29/32	LMS 61165604	LMS 61165614	10
1000 ml	NS 29/32	LMS 61165605	LMS 61165615	2
2000 ml	NS 29/32	LMS 61165606	LMS 61165616	2



SEPARATING FUNNELS

acc. to Squibb acc. to DIN ISO 4800, with plastic stopper, graduation,
borosilicate glass


Capacity	Neck size	PTFE stopcock key	glass stopcock key	Pack of
50 ml	NS 19/26	LMS 61165620	LMS 61165101	10
100 ml	NS 19/26	LMS 61165621	LMS 61165102	10
250 ml	NS 29/32	LMS 61165622	LMS 61165103	10
500 ml	NS 29/32	LMS 61165623	LMS 61165104	10
1000 ml	NS 29/32	LMS 61165624	LMS 61165105	2
2000 ml	NS 29/32	LMS 61165625	LMS 61165106	2



DROPPING FUNNELS

cylindrical form, acc. to DIN ISO 4800, with conical ground socket and cone,
plastic stopper, graduation, borosilicate glass

Capacity	NS size	PTFE stopcock key	glass stopcock key	Pack of
50 ml	NS 14/23	LMS 61164551	LMS 61164531	1
100 ml	NS 14/23	LMS 61164552	LMS 61164532	10
250 ml	NS 29/32	LMS 61164553	LMS 61164534	10
500 ml	NS 29/32	LMS 61164554	LMS 61164535	10
1000 ml	NS 29/32	LMS 61164555	LMS 61164536	10
2000 ml	NS 29/32	LMS 61164556	LMS 61164537	2



BURETTES


acc. to Mohr, with straight interchangeable glass stopcock, boro 3.3.

Capacity	Subdiv.ml	CL. AS	Pack of
10 ml	0,02	LMS 61346411	2
25 ml	0,05	LMS 61346413	2
50 ml	0,1	LMS 61346415	2
100 ml	0,2	LMS 61346416	2

SCHELLBACH BURETTES

acc. to Mohr, with straight interchangeable glass stopcock, boro 3.3.

Capacity	Subdiv.ml	CL. AS	Pack of
10 ml	0,02	LMS 61346441	2
25 ml	0,05	LMS 61346443	2
50 ml	0,1	LMS 61346445	2
100 ml	0,2	LMS 61346446	2



BURETTES


acc. to Mohr, DIN ISO 385, with straight PTFE stopcock, boro 3.3.

Capacity	Subdiv.ml	CL. AS	Pack of
10 ml	0,02	LMS 61346511	2
25 ml	0,05	LMS 61346513	2
50 ml	0,1	LMS 61346515	2
100 ml	0,2	LMS 61346516	2

SCHELLBACH BURETTES

acc. to Mohr, DIN ISO 385, with straight PTFE stopcock, boro 3.3.

Capacity	Subdiv.ml	CL. AS	Pack of
10 ml	0,02	LMS 61346541	2
25 ml	0,05	LMS 61346543	2
50 ml	0,1	LMS 61346545	2
100 ml	0,2	LMS 61346546	2



BURETTES

acc. to Mohr, DIN ISO 385 with straight PTFE spindle stopcock, boro 3.3.

Capacity	Subdiv.ml	CL. AS	Pack of
10 ml	0,02	LMS 61346611	2
25 ml	0,05	LMS 61346613	2
50 ml	0,1	LMS 61346615	2
100 ml	0,2	LMS 61346616	2

SCHELLBACH BURETTES

acc. to Mohr, DIN ISO 385, with straight PTFE spindle stopcock, boro 3.3.


Capacity	Subdiv.ml	CL. AS	Pack of
10 ml	0,02	LMS 61346641	2
25 ml	0,05	LMS 61346643	2
50 ml	0,1	LMS 61346645	2
100 ml	0,2	LMS 61346646	2



BURETTES

acc. to Daffert, Cl. AS, Schellbach, blue graduation, stopcock with two oblique bores, borosilicate glass

Capacity	Subdiv. ml	interchangeable glass key	PTFE key	Pack of
25	0,05 ml	LMS 61447401	LMS 61447400	4
50	0,1 ml	LMS 61447405	LMS 61447404	4
100	0,2 ml	LMS 61447407	LMS 61447406	4




MICROBURETTES

acc. to Bang, with straight glass stopcock Cl. B, with glass key, borosilicate glass

Capacity	Subdiv. ml	clear glass	Schellbach	Pack of
1 ml	0,01	LMS 61342211	LMS 61342261	4
2 ml	0,01	LMS 61342212	LMS 61342262	4
5 ml	0,01	LMS 61342214	LMS 61342264	4
10 ml	0,02	LMS 61342215	LMS 61342265	4

acc. to Bang, with straight glass stopcock Cl. AS, with glass key, borosilicate glass

Capacity	Subdiv. ml	clear glass	Schellbach	Pack of
1 ml	0,01	LMS 61342231	LMS 61342281	4
2 ml	0,01	LMS 61342232	LMS 61342282	4
5 ml	0,01	LMS 61342234	LMS 61342284	4
10 ml	0,02	LMS 61342235	LMS 61342285	4




MICROBURETTES

acc. to Bang, with straight PTFE spindle stopcock,
Cl. B, with PTFE intermedium stopcock, borosilicate glass

Capacity	Subdiv. ml	clear glass	Schellbach	Pack of
1 ml	0,01	LMS 61342611	LMS 61342661	4
2 ml	0,01	LMS 61342612	LMS 61342662	4
5 ml	0,01	LMS 61342614	LMS 61342664	4
10 ml	0,02	LMS 61342615	LMS 61342665	4

acc. to Bang, with straight PTFE spindle stopcock
Cl. AS, with PTFE intermedium stopcock, borosilicate glass

Capacity	Subdiv. ml	clear glass	Schellbach	Pack of
1 ml	0,01	LMS 61342631	LMS 61342481	4
2 ml	0,01	LMS 61342632	LMS 61342482	4
5 ml	0,01	LMS 61342634	LMS 61342484	4
10 ml	0,02	LMS 61342635	LMS 61342485	4




AUTOMATIC BURETTES

acc. to Dr. Schilling, Cl. B, Schellbach with automatic zero adjustment, for rapid titration with pinchcock for delivery drop by drop, blue graduation, consisting of: burette, plastic bottle and plastic base

Capacity	code	Pack of
5 ml	LMS 61046401	10
10 ml	LMS 61046402	10
25 ml	LMS 61046403	10
50 ml	LMS 61046404	10


Spare burettes

Capacity	code	Pack of
5 ml	LMS 61046411	10
10 ml	LMS 61046412	10
25 ml	LMS 61046413	10
50 ml	LMS 61046414	10



AUTOMATIC BURETTES

acc. to Pellet, without intermediate stopcock, Cl. B, with glass key, without bottle, borosilicate glass



Capacity	Subdiv. ml	clear glass	Schellbach	Pack of
10 ml	0,02	LMS 61347401	LMS 61347431	2
25 ml	0,05	LMS 61347403	LMS 61347433	2
50 ml	0,1	LMS 61347405	LMS 61347435	2


AUTOMATIC BURETTES

acc. to Pellet, without intermediate stopcock, Cl. AS, with glass key, without bottle, borosilicate glass

Capacity	Subdiv. ml	clear glass	Schellbach	Pack of
10 ml	0,02	LMS 61347411	LMS 61347441	2
25 ml	0,05	LMS 61347413	LMS 61347443	2
50 ml	0,1	LMS 61347415	LMS 61347445	2

AUTOMATIC BURETTES

acc. to Pellet, without intermediate stopcock, Cl. B, with PTFE spindle stopcock, without bottle, borosilicate glass



Capacity	Subdiv. ml	clear glass	Schellbach	Pack of
10 ml	0,02	LMS 61347601	LMS 61347631	2
25 ml	0,05	LMS 61347603	LMS 61347633	2
50 ml	0,1	LMS 61347605	LMS 61347635	2


AUTOMATIC BURETTES

acc. to Pellet, without intermediate stopcock, Cl. AS, with PTFE spindle stopcock, without bottle, borosilicate glass

Capacity	Subdiv. ml	clear glass	Schellbach	Pack of
10 ml	0,02	LMS 61347611	LMS 61347641	2
25 ml	0,05	LMS 61347613	LMS 61347643	2
50 ml	0,1	LMS 61347615	LMS 61347645	2

AUTOMATIC BURETTES

acc. to Pellet, without intermediate stopcock, Cl. B, with PTFE key, without bottle, borosilicate glass



Capacity	Subdiv. ml	clear glass	Schellbach	Pack of
10 ml	0,02	LMS 61347501	LMS 61347531	2
25 ml	0,05	LMS 61347503	LMS 61347533	2
50 ml	0,1	LMS 61347505	LMS 61347535	2

AUTOMATIC BURETTES

acc. to Pellet, without intermediate stopcock, Cl. AS, with PTFE key, without bottle, borosilicate glass

Capacity	Subdiv. ml	clear glass	Schellbach	Pack of
10 ml	0,02	LMS 61347511	LMS 61347541	2
25 ml	0,05	LMS 61347513	LMS 61347543	2
50 ml	0,1	LMS 61347515	LMS 61347545	2

RESERVOIR BOTTLES

for Automatic Burettes, NS 29/32


Capacity	clear glass	amber glass	Pack of
2000 ml	LMS 63231953	LMS 63231952	4



RED RUBBER BULBS

for Automatic Burettes


Type	Code	Pack of
Single bulb	BOE 110	1



LIEBIG CONDENSERS

borosilicate glass, hose connection diam. 11 mm


Jacket length	Cone/socket size	Code	Pack of
100 mm	NS 14/23	LMS 61451118	4
160 mm	NS 14/23	LMS 61451119	4
160 mm	NS 19/26	LMS 61451120	4
200 mm	NS 14/23	LMS 61451121	4
200 mm	NS 29/32	LMS 61451122	4
250 mm	NS 19/26	LMS 61451123	2
250 mm	NS 24/29	LMS 61451124-40	2
250 mm	NS 24/40	LMS 61451224	2
250 mm	NS 29/32	LMS 61451125	2
400 mm	NS 24/29	LMS 61451202	1
400 mm	NS 24/40	LMS 61451202-40	1
400 mm	NS 29/32	LMS 61451203	1



COILED CONDENSERS

borosilicate glass


Jacket length	Cone/socket size	Code	Pack of
200 mm	NS 14/23	LMS 61455121	4
200 mm	NS 29/32	LMS 61455122	4
300 mm	NS 29/32	LMS 61455132	2
400 mm	NS 24/40	LMS 61455930	1
400 mm	NS 29/32	LMS 61455142	1




ALLIHIN CONDENSERS

borosilicate glass, hose connection diam. 11 mm

Jacket length	Cone/socket size	Code	Pack of
160 mm	NS 14/23	LMS 61452129	4
160 mm	NS 19/26	LMS 61452130	4
250 mm	NS 19/26	LMS 61452131	2
250 mm	NS 24/29	LMS 61452132	2
250 mm	NS 24/40	LMS 61452225	2
250 mm	NS 29/32	LMS 61452133	2
300 mm	NS 29/32	LMS 61452134	1
400 mm	NS 24/40	LMS 61452134-40	1
400 mm	NS 29/32	LMS 61452135	1






DIMROTH CONDENSERS

borosilicate glass, hose connection diam. 11 mm


Jacket length	Cone/socket size	Code	Pack of
160 mm	NS 14/23	LMS 61453120	4
160 mm	NS 19/26	LMS 61453121	4
250 mm	NS 19/26	LMS 61453122	2
250 mm	NS 24/29	LMS 61453123	2
250 mm	NS 29/32	LMS 61453124	2
400 mm	NS 24/29	LMS 61453125	1
400 mm	NS 29/32	LMS 61453126	1



JACKET COILED CONDENSERS

borosilicate glass


Jacket length	Cone/socket size	Code	Pack of
200 mm	NS 29/32	LMS 61454119	2
250 mm	NS 29/32	LMS 61454120	2
300 mm	NS 29/32	LMS 61454121	2
400 mm	NS 29/32	LMS 61454122	2



SOXHLET EXTRACTORS

cones NS 29/32, (NS 24/40 for L version) borosilicate glass

Vol. extractor	socket size	Code	Pack of
30 ml	NS 29/32	LMS 61435100	10
50 ml	NS 29/32	LMS 61435101	10
100 ml	NS 45/40	LMS 61435102	10
100 ml	NS 24/40	LMS 61435102L	10
150 ml	NS 45/40	LMS 61435103	10
250 ml	NS 45/40	LMS 61435104	10
250 ml	NS 24/40	LMS 61435104L	10
500 ml	NS 60/46	LMS 61435105	10




SOXHLET EXTRACTION APPARATUS

complete with: Soxhlet extractor, flat bottom flask NS 29/32, Dimroth condenser NS 29/32, reduction adapter

Vol. extractor	flat bottom flask	Code	Pack of
30 ml	100 ml	LMS 61612102	1
50 ml	250 ml	LMS 61612103	1
100 ml	250 ml	LMS 61612104	1
150 ml	500 ml	LMS 61612105	1
250 ml	500 ml	LMS 61612106	1
500 ml	1000 ml	LMS 61612107	1
500 ml	2000 ml	LMS 616121072	1

WEIGHING BOTTLES

with interchangeable knob lid, matt label, low form



Ø x Height	Code	Pack of
25 x 25 mm	LMS 61135102	10
35 x 30 mm	LMS 61135106	10
45 x 30 mm	LMS 61135109	10
50 x 30 mm	LMS 61135111	10
60 x 30 mm	LMS 61135113	10
80 x 30 mm	LMS 61135116	10

WEIGHING BOTTLES

with interchangeable knob lid, matt label, tall form


Ø x Height	Code	Pack of
25 x 40 mm	LMS 61135104	10
30 x 50 mm	LMS 61135107	10
35 x 70 mm	LMS 61135137	10
50 x 80 mm	LMS 61135158	10
60 x 80 mm	LMS 61135168	10



ERLENMEYER FLASKS

with conical ground socket, borosilicate glass 3.3., with approximated volume scale, acc. to DIN 4797

Capacity	Socket size	Code	Pack of
25 ml	NS 14 / 23	LMS 61217601	10
25 ml	NS 19 / 26	LMS 61217604	10
50 ml	NS 14 / 23	LMS 61217602	10
50 ml	NS 19 / 26	LMS 61217605	10
50 ml	NS 24 / 29	LMS 61217611	10
50 ml	NS 29 / 32	LMS 61217621	10
100 ml	NS 14 / 23	LMS 61217603	10
100 ml	NS 19 / 26	LMS 61217606	10
100 ml	NS 24 / 29	LMS 61217612	10
100 ml	NS 29 / 32	LMS 61217622	10
200 ml	NS 19 / 26	LMS 61217607	10
200 ml	NS 24 / 29	LMS 61217613	10
200 ml	NS 29 / 32	LMS 61217623	10
200 ml	NS 45 / 40	LMS 61217629	10
250 ml	NS 19 / 26	LMS 61217608	10
250 ml	NS 24 / 29	LMS 61217614	10
250 ml	NS 29 / 32	LMS 61217624	10
250 ml	NS 45 / 40	LMS 61217630	10
300 ml	NS 19 / 26	LMS 61217609	10
300 ml	NS 24 / 29	LMS 61217615	10
300 ml	NS 29 / 32	LMS 61217625	10
300 ml	NS 45 / 40	LMS 61217631	10
500 ml	NS 24 / 29	LMS 61217616	10
500 ml	NS 29 / 32	LMS 61217626	10
500 ml	NS 45 / 40	LMS 61217632	10
1000 ml	NS 24 / 29	LMS 61217617	10
1000 ml	NS 29 / 32	LMS 61217627	10
1000 ml	NS 45 / 40	LMS 61217633	10
2000 ml	NS 24 / 29	LMS 61217618	10
2000 ml	NS 29 / 32	LMS 61217628	10
2000 ml	NS 45 / 40	LMS 61217634	10



STOPPERS

with flat and pointed bottom, borosilicate glass

Size	Code	Pack of
NS 5 / 13	LMS 61424401	25
NS 7 / 16	LMS 61424402	25
NS 10 / 19	LMS 61424403	25
NS 12 / 21	LMS 61424404	25
NS 14 / 23	LMS 61424405	25
NS 19 / 26	LMS 61424406	25
NS 24 / 29	LMS 61424407	10
NS 29 / 32	LMS 61424408	10



FLAT BOTTOM FLASKS

with conical ground socket, borosilicate glass 3.3.



Capacity	Socket size	Code	Pack of
25 ml	NS 14 / 23	LMS 61212201	10
50 ml	NS 14 / 23	LMS 61212202	10
50 ml	NS 19 / 26	LMS 61212204	10
50 ml	NS 24 / 29	LMS 61212207	10
50 ml	NS 29 / 32	LMS 61212215	10
100 ml	NS 14 / 23	LMS 61212203	10
100 ml	NS 19 / 26	LMS 61212205	10
100 ml	NS 24 / 29	LMS 61212208	10
100 ml	NS 29 / 32	LMS 61212216	10
250 ml	NS 19 / 26	LMS 61212206	10
250 ml	NS 24 / 29	LMS 61212209	10
250 ml	NS 24 / 40	LMS 61212212-40	10
250 ml	NS 29 / 32	LMS 61212217	10
500 ml	NS 24 / 29	LMS 61212210	10
500 ml	NS 24 / 40	LMS 61212213-40	10
500 ml	NS 29 / 32	LMS 61212218	10
1000 ml	NS 24 / 29	LMS 61212211	10
1000 ml	NS 24 / 40	LMS 61212214-40	10
1000 ml	NS 29 / 32	LMS 61212219	10
1000 ml	NS 45 / 40	LMS 61212221	10
2000 ml	NS 45 / 40	LMS 61212222	10

ROUND BOTTOM FLASKS

with conical ground socket, borosilicate glass 3.3.



Capacity	Socket size	Code	Pack of
10 ml	NS 14 / 23	LMS 61211200	10
25 ml	NS 14 / 23	LMS 61211201	10
25 ml	NS 19 / 26	LMS 61211205	10
50 ml	NS 14 / 23	LMS 61211202	10
50 ml	NS 19 / 26	LMS 61211206	10
50 ml	NS 24 / 29	LMS 61211210	10
50 ml	NS 29 / 32	LMS 61211216	10
100 ml	NS 14 / 23	LMS 61211203	10
100 ml	NS 19 / 26	LMS 61211207	10
100 ml	NS 24 / 29	LMS 61211211	10
100 ml	NS 29 / 32	LMS 61211217	10
250 ml	NS 14 / 23	LMS 61211204	10
250 ml	NS 19 / 26	LMS 61211208	10
250 ml	NS 24 / 29	LMS 61211212	10
250 ml	NS 24 / 40	LMS 61211212-40	10
250 ml	NS 29 / 32	LMS 61211218	10
250 ml	NS 34 / 35	LMS 61211228	10
250 ml	NS 45 / 40	LMS 61211236	10
500 ml	NS 19 / 26	LMS 61211209	10
500 ml	NS 24 / 29	LMS 61211213	10
500 ml	NS 24 / 40	LMS 61211213-40	10
500 ml	NS 29 / 32	LMS 61211219	10
500 ml	NS 34 / 35	LMS 61211229	10
500 ml	NS 45 / 40	LMS 61211237	10
1000 ml	NS 24 / 29	LMS 61211214	10
1000 ml	NS 24 / 40	LMS 61211214-40	10
1000 ml	NS 29 / 32	LMS 61211220	10
1000 ml	NS 34 / 35	LMS 61211230	10
1000 ml	NS 45 / 40	LMS 61211238	10
2000 ml	NS 24 / 29	LMS 61211215	10
2000 ml	NS 29 / 32	LMS 61211221	10
2000 ml	NS 34 / 35	LMS 61211231	10
2000 ml	NS 45 / 40	LMS 61211239	10

JIPO PORCELAIN

CRUCIBLES, LOW FORM

Ø (mm)	H (mm)	Vol (ml)	Code	Pack of
30	19	5	JIP 100301	5
35	22	10	JIP 100351	5
40	25	17	JIP 100401	5
45	28	21	JIP 100451	5
50	32	34	JIP 100501	5
60	38	62	JIP 100601	5
70	44	91	JIP 100701	5



CRUCIBLES, MIDDLE FORM

Ø (mm)	H (mm)	Vol (ml)	Code	Pack of
30	25	10	JIP 200301	5
35	28	12	JIP 200351	5
40	32	20	JIP 200401	5
45	36	30	JIP 200451	5
50	40	45	JIP 200501	5
60	48	80	JIP 200601	5
70	56	120	JIP 200701	5
80	82	200	JIP 200801	5



CRUCIBLES, HIGH FORM

Ø (mm)	H (mm)	Vol (ml)	Code	Pack of
30	38	15	JIP 300301	5
35	44	26	JIP 300351	5
40	50	35	JIP 300401	5
45	56	50	JIP 300451	5
50	62	72	JIP 300501	5
60	75	120	JIP 300601	5



LIDS FOR CRUCIBLES

Ø (mm) crucible	Ø (mm) lid	Code	Pack of
30	34	JIP 000301	5
35	39	JIP 000351	5
40	44	JIP 000401	5
45	49	JIP 000451	5
50	54	JIP 000501	5
60	64	JIP 000601	5
70	74	JIP 000701	5
80	84	JIP 000801	5



GOOCH CRUCIBLES

Ø (mm)	H (mm)	Vol (ml)	Code	Pack of
30	36	17	JIP 234101	5
35	43	25	JIP 234301	5
39	49	35	JIP 234401	5
60	71	120	JIP 234501	5



FILTER CRUCIBLES

with porous base



H (mm)	Vol (ml)	Porosity	Code	Pack of
28	8	P1 = 6μ	JIP 400256	2
28	8	P2 = 7μ	JIP 400257	2
28	8	P3 = 8μ	JIP 400258	2
35	15	P1 = 6μ	JIP 400306	2
35	15	P2 = 7μ	JIP 400307	2
35	15	P3 = 8μ	JIP 400308	2
40	25	P1 = 6μ	JIP 400356	2
40	25	P2 = 7μ	JIP 400357	2
40	25	P3 = 8μ	JIP 400358	2
43	35	P1 = 6μ	JIP 400406	2
43	35	P2 = 7μ	JIP 400407	2
43	35	P3 = 8μ	JIP 400408	2
50	50	P1 = 6μ	JIP 400456	2
50	50	P2 = 7μ	JIP 400457	2
50	50	P3 = 8μ	JIP 400458	2

BUECHNER FUNNELS



Nom. size	Stem length	Hole Ø	Code	Pack of
45	43	1	JIP 237101	1
55	64	1	JIP 237201	1
70	64	1	JIP 237301	1
90	71	2	JIP 237401	1
110	83	2	JIP 237501	1
125	85	2	JIP 237601	1
150	96	2	JIP 237701	1
185	117	2	JIP 237801	1
240	117	2	JIP 237901	1
270	130	2	JIP 237921	1
320	133	2	JIP 237991	1

DESICCATOR PLATES, ACC. TO CSN



Ø (mm)	for Desic.	Code	Pack of
90	100	JIP 238101	1
140	150	JIP 238331	1
189	200	JIP 238521	1
241	250	JIP 238821	1
290	300	JIP 238991	1

DESICCATOR PLATES, ACC. TO DIN



Ø (mm)	for Desic.	Code	Pack of
90	100	JIP 138101	1
140	150	JIP 138201	1
190	200	JIP 138301	1
235	250	JIP 138401	1
280	300	JIP 138501	1

EVAPORATION DISH, ROUND BOTTOM

deep form, with lip

Ø (mm)	H (mm)	Vol (ml)	Code	Pack of
40	16	10	JIP 206001	5
50	20	20	JIP 206101	5
63	25	30	JIP 206121	5
80	32	60	JIP 206201	5
81	38	107	JIP 206221	5
94	42	140	JIP 206301	5
100	40	150	JIP 206321	5
112	47	232	JIP 206401	5
125	50	285	JIP 206411	5
132	55	300	JIP 206421	5
138	61	400	JIP 206501	5
160	64	580	JIP 206611	5
170	74	700	JIP 206621	1
200	80	1000	JIP 206721	1
202	86	1300	JIP 206801	1
320	140	6000	JIP 206991	1



EVAPORATION DISH, ROUND BOTTOM

semideep form, with lip

Ø (mm)	H (mm)	Vol (ml)	Code	Pack of
72	26	41	JIP 207101	5
83	27	62	JIP 207201	5
96	30	110	JIP 207301	5
110	37	160	JIP 207401	5
117	41	215	JIP 207421	5
127	42	240	JIP 207501	5
142	46	325	JIP 207521	5
154	50	420	JIP 207601	5
170	55	558	JIP 207701	5



EVAPORATION DISH, FLAT BOTTOM

with lip

Ø (mm)	H (mm)	Vol (ml)	Code	Pack of
54	22	22	JIP 274101	5
70	30	62	JIP 274121	5
86	33	93	JIP 274201	5
98	40	154	JIP 274301	5
112	50	265	JIP 274401	5
126	53	330	JIP 274501	5
145	55	450	JIP 274601	5
151	63	620	JIP 274701	5
169	68	800	JIP 274801	1
195	74	1000	JIP 274901	1
254	105	2500	JIP 274991	1



MORTARS WITH LIP, GLAZED



Ø (mm)	H (mm)	Vol (ml)	Code	Pack of
51	32	25	JIP 211001	5
56	36	30	JIP 211101	2
63	41	70	JIP 211111	2
76	45	75	JIP 211121	2
83	50	110	JIP 211201	2
90	56	160	JIP 211221	1
105	64	220	JIP 211301	1
125	71	400	JIP 211321	1
139	82	500	JIP 211401	1
150	90	650	JIP 211501	1
180	104	1000	JIP 211601	1
224	127	2250	JIP 211701	1
268	150	4300	JIP 211801	1

PESTLES, GLAZED



Ø (mm)	length (mm)	Code	Pack of
14	54	JIP 213001	5
24	115	JIP 213101	2
30	135	JIP 213201	2
36	150	JIP 213301	2
42	175	JIP 213401	1
55	210	JIP 213501	1
74	250	JIP 213601	1
85	280	JIP 213701	1
101	370	JIP 213801	1

MORTARS WITH LIP, ROUGH



Ø (mm)	H (mm)	Vol (ml)	Code	Pack of
51	32	25	JIP 211000	5
56	36	30	JIP 211100	2
63	41	70	JIP 211110	2
76	45	75	JIP 211120	2
83	50	110	JIP 211200	2
90	56	160	JIP 211220	1
105	64	220	JIP 211300	1
125	71	400	JIP 211320	1
139	82	500	JIP 211400	1
150	90	650	JIP 211500	1
180	104	1000	JIP 211600	1
224	127	2250	JIP 211700	1
268	160	4300	JIP 211800	1

PESTLES, ROUGH



Ø (mm)	length (mm)	Code	Pack of
14	54	JIP 213000	5
24	115	JIP 213100	2
30	135	JIP 213200	2
36	150	JIP 213300	2
42	175	JIP 213400	1
55	210	JIP 213500	1
74	250	JIP 213600	1
85	280	JIP 213700	1
101	370	JIP 213800	1

BOECO

The family-owned company Boeckel + Co. (GmbH + Co.) KG was founded in 1929 and now exists in the fourth generation.

To our customers on all 5 continents we are pleased to introduce our fifth general BOECO catalog with our extended and updated complete line of laboratory equipment and consumables.

We offer high quality at competitive prices by a convenient and cost-effective one-stop supply.

CONTACT

Boeckel + Co. (GmbH + Co.) KG
Roedingsmarkt 33
20459 Hamburg
Germany

Phone: + 49 40 325627-0
Fax: + 49 40 325627-41
E-mail: info@boeco.com
Internet: www.boeco.com

V.A.T. ID-No: DE 118903538

BANK

EURO €

Deutsche Bank AG, Hamburg,
SWIFT: DEUTDEHHXXX
IBAN- No. DE 65 2007 0000 0010 5601 00

USD \$

Deutsche Bank AG, Hamburg,
SWIFT: DEUTDEHHXXX
IBAN- No. DE 27 2007 0000 0010 5601 05