


PIPETTE TIPS

Cat. No.	Volume	Length	Type	Description	Packaging
					
BOE 1-300-81-0	0,1-10 µl	31,2 mm	A	Top Line Natural Tip	bag / 1000 pcs. / carton 20.000 pcs.
BOE 1-300-96-0	0,1-10 µl		A	Top Line Natural Tip	96-tip Rackbox / pack 5x10 racks
BOE 1-301-96-0	0,1-10 µl		A	Top Line Natural Tip, sterile	96-tip Rackbox / pack 5x10 racks
BOE 2-301-96-0	0,1-10 µl		A	Top Line Filter Tip, sterile	96-tip Rackbox / pack 5x10 racks
					
BOE 03.033.0100	0,1-20 µl	40,5 mm	B	Standard Tip, Natural, long type	bag / 1000 pcs. / carton / 10.000 pcs.
					
BOE 1-350-81-0	0,5-10/20 µl	45,7 mm	C	Top Line Natural Tip, long type	bag / 1000 pcs. / carton / 20.000 pcs.
BOE 1-350-96-0	0,5-10/20 µl		C	Top Line Natural Tip, long type	96-tip Rackbox / pack 5x10 racks
BOE 1-351-96-0	0,5-10/20 µl		C	Top Line Natural Tip, long type, sterile	96-tip Rackbox / pack 5x10 racks
					
BOE 03.013.0103	2- 200 µl	52,6 mm	D	Standard Tip, yellow	bag / 1000 pcs. / carton / 25.000 pcs.
 					
BOE 1-410-81-1	1-200 µl	53,4 mm	E	Top Line Yellow Tip, bevelled, graduated	bag / 1000 pcs. / carton / 30.000 pcs.
BOE 1-410-96-1	1-200 µl		E	Top Line Yellow Tip, bevelled graduated, nonsterile	96-tip Rackbox / pack 5x10 racks
BOE 1-411-96-0	1-200 µl		E	Top Line Natural Tip, bevelled graduated, sterile	96-tip Rackbox / pack 5x10 racks
BOE 2-416-96-0	1-20/30 µl		E	Top Line Filter Tip, natural, sterile	96-tip Rackbox / pack 5x10 racks
BOE 2-412-96-0	1-100 µl		E	Top Line Filter Tip, natural, sterile	96-tip Rackbox / pack 5x10 racks
BOE 2-413-96-0	1-200 µl		E	Top Line Filter Tip, natural, sterile	96-tip Rackbox / pack 5x10 racks
BOE 2-415-96-0	5-300 µl		E	Top Line Filter Tip, natural, sterile	96-tip Rackbox / pack 5x10 racks
					
BOE 732031	5-300 µl	52 mm	F	Natural Tip	box / 1000 pcs. / carton / 10.000 pcs.
					
BOE 03.023.0101	50-1000 µl	70,5 mm	G	Standard Tip, blue	bag / 1000 pcs. / carton / 10.000 pcs.
 					
BOE 1-600-80-2	100-1000 µl	84,3 mm	I	Top Line Tip, blue, bevelled, graduated	bag / 500 pcs. / carton / 10.000 pcs.
BOE 1-601-96-2	100-1000 µl		I	Top Line Tip, blue, bevelled graduated	96-tip Rackbox / pack 4 x 8 racks
BOE 1-601-96-0	100-1000 µl		I	Top Line Tip, clear, bevelled graduated, sterile	96-tip Rackbox / pack 4 x 8 racks
BOE 2-601-96-0	100-1000 µl		I	Top Line Filter Tip, sterile	96-tip Rackbox / pack 4 x 8 racks

BOECO Top Line Tips are certified to be free of any detectable human DNA, DNase, RNase and PCR-inhibitors

Cat.No.	Volume	Length	Type	Description	Packaging
---------	--------	--------	------	-------------	-----------



BOE 9600650	1 - 5 ml	126 mm	J	Natural Tip, for SA/IP Series	carton / 300 pcs.
-------------	----------	--------	---	-------------------------------	-------------------

This Tip is also suitable for our previous Standard and SP pipettes with 2000 µl and 5000 µl as well as BIOHIT® pipettes



BOE 702601	0,5 - 5 ml	160 mm	K	Natural Tip, for GP Series	bag / 200 pcs. / carton 1.000 pcs.
------------	------------	--------	---	----------------------------	------------------------------------

This Tip is also suitable for Thermo Fisher Scientific Finn timer® F1 and F2 of 5 ml



BOE 1218610	1 - 5 ml	124 mm	N	Natural Tip, for IP Series	bag / 250 pcs. / carton 2.500 pcs.
-------------	----------	--------	---	----------------------------	------------------------------------

This Tip is also suitable for Gilson® and HTL® pipettes of 5 ml



BOE 9600640	2 - 10 ml	150 mm	L	Natural Tip, for SA/IP Series	bag / 100 pcs.
-------------	-----------	--------	---	-------------------------------	----------------

This Tip is also suitable for Thermo Fisher Scientific Finn timer® F1 and F2 of 10 ml



BOE 702604	1 - 10 ml	156,5 mm	M	Natural Tip, for GP Series	bag / 100 pcs. / carton 1.000 pcs.
------------	-----------	----------	---	----------------------------	------------------------------------

This Tip is also suitable for the old SP Pipettes with 10 ml volume, as well as for Eppendorf® and Gilson® pipettes of 10 ml



BOE 1220500	1 - 10 ml	181 mm	O	Natural Tip, for IP Series	bag / 200 pcs. / carton 1.000 pcs.
-------------	-----------	--------	---	----------------------------	------------------------------------

This Tip is also suitable for Gilson® and Eppendorf® pipettes of 10 ml



BOECO RACKBOXES

BOECO

The family-owned company
Boeckel + Co. (GmbH + Co.) KG was founded in
1929 and now exists in the fourth generation.

To our customers on all 5 continents we are
pleased to introduce our fifth general BOECO
catalog with our extended and updated
complete line of laboratory equipment and
consumables.

We offer high quality at competitive prices by
a convenient and cost-effective one-stop supply.

CONTACT

Boeckel + Co. (GmbH + Co.) KG
Roedingsmarkt 33
20459 Hamburg
Germany

Phone: + 49 40 325627-0
Fax: + 49 40 325627-41
E-mail: info@boeco.com
Internet: www.boeco.com

V.A.T. ID-No: DE 118903538

BANK

EURO €

Deutsche Bank AG, Hamburg,
SWIFT: DEUTDEHHXXX
IBAN- No. DE 65 2007 0000 0010 5601 00

USD \$

Deutsche Bank AG, Hamburg,
SWIFT: DEUTDEHHXXX
IBAN- No. DE 27 2007 0000 0010 5601 05