



**NARROW NECK BOTTLE
WITH SPRAYING TUBE**

NARROW-NECK / WASH BOTTLES

Narrow neck bottles in polyethylene soft (LD-PE), graduated, natural

Volume	Code	Dimensions D-d-H (mm)	Pack
50 ml	KJS E50	37 - 11 - 87	100
100 ml	KJS E100	47 - 11 - 105	550
250 ml	KJS E250	63 - 18 - 119	252
500 ml	KJS E500	78 - 18 - 152	123
1000 ml	KJS E1000	96 - 20 - 206	61
2000 ml	KJS E2000	120 - 20 - 264	36

Screw caps for wash bottles, with spraying tube (swan-neck)

Diam.	Code	fit to narrow neck bottles:	Pack
18 mm	KJS SR18	E50, E100	25
25 mm	KJS SR25	E250, E500	25
28 mm	KJS SR28	E1000, E2000	25



WIDE-NECK BOTTLES

round, in polyethylene soft (LD-PE), natural

Volume	Code	Dimensions D-d-H (mm)	Pack
50 ml	KJS W50	39 - 24 - 77	100
100 ml	KJS W100	47 - 24 - 94	624
250 ml	KJS W250	63 - 31 - 119	279
500 ml	KJS W500	78 - 39 - 152	141
1000 ml	KJS W1000	96 - 50 - 206	64
1500 ml	KJS W1500	109 - 50 - 226	45
2000 ml	KJS W2000	120 - 50 - 247	33

Screw caps with threads for narrow- and wide-neck bottles (LD-PE)

Diam.	Code	fit to bottles:	Pack
18 mm	KJS V18	E50, E100	200
25 mm	KJS V25	E250, E500	200
28 mm	KJS V28	E1000, E2000	200
32 mm	KJS V32	W50, W100	200
40 mm	KJS V40	W250	100
50 mm	KJS V50	W500	100
65 mm	KJS V65	W1000, W1500, W2000	100



SQUARE BOTTLES, WIDE NECK

made in HD-PE, natural, with black tamper evident cap packed aside

Volume	Code	Dimensions (mm)	Pack
100 ml	KJS 141-0100	52 x 52 mm, height 82 mm	440
250 ml	KJS 141-0250	67 x 67 mm, height 107 mm	215
500 ml	KJS 141-0250	85 x 85 mm, height 129 mm	123
1000 ml	KJS 141-01000	85 x 85 mm, height 210 mm	77



TEST TUBE RACKS

for general purpose, square holes, easy to assemble,
white, operating temperature -20 to +90 °C

Base dimension: 245 x 105 mm


Tube Ø	Code	Capacity
13	BOE 290131	90
16	BOE 260161	60
20	BOE 240201	40
25	BOE 224251	24
30	BOE 221301	21

PIPETTE RACK, PP

Rotary model for 94 pipettes of different sizes,
autoclavable at temperatures up to 121°C.

Code

BOE 194



MAGNETIC PTFE ENCAPSULATED STIRRING BARS

cylindrical, smooth round profile, universal application

ØxLmm	Code
4,5 x 12	BOE 114512
6 x 10	BOE 110610
6 x 13	BOE 110613
6 x 20	BOE 110620
6 x 25	BOE 110625
7 x 20	BOE 110720
7 x 30	BOE 110730
7 x 40	BOE 110740
8 x 50	BOE 110850
10 x 60	BOE 111060
10 x 70	BOE 111070
10 x 80	BOE 111080

MAGNETIC PTFE ENCAPSULATED STIRRING BARS

with pivot ring

ØxLmm	Code
6 x 25	BOE 130625
8 x 40	BOE 130840
8 x 50	BOE 130850
10 x 25	BOE 131025
10 x 60	BOE 131060

MAGNETIC STIRRING BAR RETRIEVER, PTFE

Length mm	Code
350	BOE 114350

