

LABORATORY TIMER



MECHANICAL TIMER

0-60 min, sec. alarm, housing colour white/blue

Code

BOE 124



BOECO DIGITAL STOPWATCH

waterproof,
Stopwatch precision 1/100 seconds first 30 min, 1 sec. up to max 24 h,
2 splits memory
Chronograph up to 23 hours, 59 minutes, 59,99 seconds
Calendar shows day, week, month and year
Daily alarm, hourly chime and snooze function
Time shows second, minute, hour, AM/PM & 12/24 switch
With batteries (button cells, 2pcs)

Code

BOE 260



DIGITAL TIMER

Count down 99 min., 59 sec., alarm when countdown reaches
"zero", magnet, silver, with battery (1 pc. button cell SR 44 1,5 V)

Code

BOE 200



BOECO DIGITAL TIMER

Count down/up timer, 99 min. 59 sec. alarm when countdown reaches
"zero", magnet, blue, with battery (AAA, 1 pc.)

Code

BOE 210

BOECO DIGITAL TIMER, JUMBO DISPLAY

82 x 68 x 20,5 mm
count up/ count down timer 23 hours, 59 min., 59 sec.
Clock (12h format)
Reminder: sound and flashing light
Casing material: ABS plastic
Fixed Support: Magnet + stand + hanging role
Battery: 2 pcs. AAA



Code

BOE 270

BOECO 2-LINE DIGITAL TIMER, JUMBO DISPLAY

92 x 79 x 21 mm
Display size: 60 x 40 mm
count up/ count down timer 99 hours, 59 min., 59 sec.
2 channels timer, stopwatch function with flashing light
Casing material: ABS plastic
Fixed support: magnet + stand + hanging hole
Battery: 2 pcs. AAA



Code

BOE 280

BOECO 3-LINE DIGITAL TIMER, JUMBO DISPLAY

3 channels digital timer
88,7 x 68 x 16,2 mm
Alarm and flashing when countdown reaches „zero“
Clock (12h format)
Clip for hanging and standing,
magnet with battery (2 x AAA, 1,5 V)



Code

BOE 290