



LABORATORY FILTER PAPER

BOECO FILTER PAPER

QUANTITATIVE FILTERS

quantitative ashless filter papers are recommended for quantitative analysis. As quantitative analysis requires a high purity, these filters are being washed out with acid



388 - black spot - fast filtering, wide pores, ashless hardened

For coarse and voluminous precipitates such as iron-aluminium and chromium hydroxide, Si-determination in steel and pig iron analysis

Basis weight: 84 g/m²
Filtration velocity: 10 ml / 10 seconds
Typical pore size: 12 - 15 µm
Similar to: S&S 589/1, M-N 640w

Diam.	Code	Grade/Description	Pack of
55 mm	FTR 3.101.055	388 / Filter discs. Black dot	100
70 mm	FTR 3.101.070	388 / Filter discs. Black dot	100
90 mm	FTR 3.101.090	388 / Filter discs. Black dot	100
110 mm	FTR 3.101.110	388 / Filter discs. Black dot	100
125 mm	FTR 3.101.125	388 / Filter discs. Black dot	100
150 mm	FTR 3.101.150	388 / Filter discs. Black dot	100
185 mm	FTR 3.101.185	388 / Filter discs. Black dot	100
240 mm	FTR 3.101.240	388 / Filter discs. Black dot	100

389 - white spot - medium fast filtering, medium wide pores, ashless hardened

Typical grade for quantitative tasks, coarser precipitates, such as lead-, iron and silver sulphide alkali carbonates, food stuff and soil analysis

Basis weight: 84 g/m²
Filtration velocity: 10 ml / 20 seconds
Typical pore size: 8 - 12 µm
Similar to: Whatman 41, S&S 589/2, M-N 640m

Diam.	Code	Grade/Description	Pack of
55 mm	FTR 3.102.055	389 / Filter discs. White dot	100
70 mm	FTR 3.102.070	389 / Filter discs. White dot	100
90 mm	FTR 3.102.090	389 / Filter discs. White dot	100
110 mm	FTR 3.102.110	389 / Filter discs. White dot	100
125 mm	FTR 3.102.125	389 / Filter discs. White dot	100
150 mm	FTR 3.102.150	389 / Filter discs. White dot	100
185 mm	FTR 3.102.185	389 / Filter discs. White dot	100
240 mm	FTR 3.102.240	389 / Filter discs. White dot	100



392 - red spot - medium to slowly filtering, medium tight, ashless hardened

Fast filtration of fine precipitates such as calcium oxalate, magnesium ammonium phosphate, coarser forms of barium sulphate

Basis weight: 84 g/m²
Filtration velocity: 10 ml / 50 seconds
Typical pore size: 5 - 8 µm
Similar to: Whatman 43, S&S 589/5, M-N 640md

Diam.	Code	Grade/Description	Pack of
55 mm	FTR 3.105.055	392 / Filter discs. Red dot	100
70 mm	FTR 3.105.070	392 / Filter discs. Red dot	100
90 mm	FTR 3.105.090	392 / Filter discs. Red dot	100
110 mm	FTR 3.105.110	392 / Filter discs. Red dot	100
125 mm	FTR 3.105.125	392 / Filter discs. Red dot	100
150 mm	FTR 3.105.150	392 / Filter discs. Red dot	100
185 mm	FTR 3.105.185	392 / Filter discs. Red dot	100
240 mm	FTR 3.105.240	392 / Filter discs. Red dot	100

390 - green spot - slowly filtering, narrow pores, ashless hardened

Filtration of fine precipitates such as lead dioxide, calcium fluoride nickel sulphide, zinc sulphide

Basis weight: 84 g/m²
Filtration velocity: 10 ml / 100 seconds
Typical pore size: 3 - 5 µm
Similar to: Whatman 40, S&S 589/6, M-N 640d

Diam.	Code	Grade/Description	Pack of
55 mm	FTR 3.103.055	390 / Filter discs. Green dot	100
70 mm	FTR 3.103.070	390 / Filter discs. Green dot	100
90 mm	FTR 3.103.090	390 / Filter discs. Green dot	100
110 mm	FTR 3.103.110	390 / Filter discs. Green dot	100
125 mm	FTR 3.103.125	390 / Filter discs. Green dot	100
150 mm	FTR 3.103.150	390 / Filter discs. Green dot	100
185 mm	FTR 3.103.185	390 / Filter discs. Green dot	100
240 mm	FTR 3.103.240	390 / Filter discs. Green dot	100



391 - blue spot - very slowly filtering, fine pores, ashless hardened

Fine-grained precipitates such as barium sulphate, metannic acid, cupros oxide

Basis weight: 84 g/m²
Filtration velocity: 10 ml / 180 seconds
Typical pore size: 2 - 3 µm
Similar to: Whatman 44, S&S 589/3, M-N 640de

Diam.	Code	Grade/Description	Pack of
55 mm	FTR 3.104.055	391 / Filter discs. Blue dot	100
70 mm	FTR 3.104.070	391 / Filter discs. Blue dot	100
90 mm	FTR 3.104.090	391 / Filter discs. Blue dot	100
110 mm	FTR 3.104.110	391 / Filter discs. Blue dot	100
125 mm	FTR 3.104.125	391 / Filter discs. Blue dot	100
150 mm	FTR 3.104.150	391 / Filter discs. Blue dot	100
185 mm	FTR 3.104.185	391 / Filter discs. Blue dot	100
240 mm	FTR 3.104.240	391 / Filter discs. Blue dot	100



393 - purple spot - particularly slow filtering, especially fine pores, ashless hardened

For extremely difficult filtration conditions and particularly fine precipitates

Basis weight: 100 g/m²
Filtration velocity: 10 ml / 250 seconds
Typical pore size: 1 - 2 µm
Similar to: Whatman 42

Diam.	Code	Grade/Description	Pack of
55 mm	FTR 3.127.055	393 / Filter discs. Purple dot	100
70 mm	FTR 3.127.070	393 / Filter discs. Purple dot	100
90 mm	FTR 3.127.090	393 / Filter discs. Purple dot	100
110 mm	FTR 3.127.110	393 / Filter discs. Purple dot	100
125 mm	FTR 3.127.125	393 / Filter discs. Purple dot	100
150 mm	FTR 3.127.150	393 / Filter discs. Purple dot	100
185 mm	FTR 3.127.185	393 / Filter discs. Purple dot	100
240 mm	FTR 3.127.240	393 / Filter discs. Purple dot	100

QUALITATIVE FILTERS

Grade 1288 - wide-pored, spongy structure, soft, very high rate of filtration

coarse-flocculent and bulky precipitates, in particular hydroxides of iron, aluminium and chromium, sulphides of copper, bismuth, cobalt and iron, various organic metal precipitates, determination of silicon in the analysis of steel and pig iron, etc.

Basis weight 84 g/m²
Filtration velocity: 10 ml / 10 seconds
Typical pore size: 12 - 15 µm
Similar to: Whatman 4, S&S 604, M-N 617

Diam.	Code	Grade/Description	Pack of
55 mm	FTR 3.206.055	1288 / Filter discs.	100
70 mm	FTR 3.206.070	1288 / Filter discs.	100
90 mm	FTR 3.206.090	1288 / Filter discs.	100
110 mm	FTR 3.206.110	1288 / Filter discs.	100
125 mm	FTR 3.206.125	1288 / Filter discs.	100
150 mm	FTR 3.206.150	1288 / Filter discs.	100
185 mm	FTR 3.206.185	1288 / Filter discs.	100
240 mm	FTR 3.206.240	1288 / Filter discs.	100

Grade 1289 - plain, medium wide pores, medium fast

commonly used grade for analytical work, coarser precipitates like sulphides of silver, arsenic, antimony, cadmium, lead, iron and magnes, lead chromate, ammonium molybdate, alkaline earth carbonates etc.

Basis weight 84 g/m²
Filtration velocity: 10 ml / 20 seconds
Typical pore size: 8 - 12 µm
Similar to: S&S 597, M-N 616

Diam.	Code	Grade/Description	Pack of
55 mm	FTR 3.207.055	1289 / Filter discs.	100
70 mm	FTR 3.207.070	1289 / Filter discs.	100
90 mm	FTR 3.207.090	1289 / Filter discs.	100
110 mm	FTR 3.207.110	1289 / Filter discs.	100
125 mm	FTR 3.207.125	1289 / Filter discs.	100
150 mm	FTR 3.207.150	1289 / Filter discs.	100
185 mm	FTR 3.207.185	1289 / Filter discs.	100
240 mm	FTR 3.207.240	1289 / Filter discs.	100

Grade 1291 - fine pored, very dense, very low rate of filtration

fine-grained precipitates cold-precipitated barium sulphate, meta-stannic acid

Basis weight 84 g/m²
Filtration velocity: 10 ml / 180 seconds
Typical pore size: 2 - 3 µm
Similar to: S&S 602 h, M-N 619 de

Diam.	Code	Grade/Description	Pack of
55 mm	FTR 3.209.055	1291 / Filter discs.	100
70 mm	FTR 3.209.070	1291 / Filter discs.	100
90 mm	FTR 3.209.090	1291 / Filter discs.	100
110 mm	FTR 3.209.110	1291 / Filter discs.	100
125 mm	FTR 3.209.125	1291 / Filter discs.	100
150 mm	FTR 3.209.150	1291 / Filter discs.	100
185 mm	FTR 3.209.185	1291 / Filter discs.	100
240 mm	FTR 3.209.240	1291 / Filter discs.	100

Grade 292 - medium to slowly filtering, medium tight

fast filtration of fine precipitates such as calcium oxalate, magnesium ammonium phosphate, coarser forms of barium sulphate

Basis weight 87 g/m²
Filtration velocity: 10 ml / 45 seconds
Typical pore size: 5 - 8 µm
Similar to: Whatman 1, S&S 593, M-N 616 md



Diam.	Code	Grade/Description	Pack of
55 mm	FTR 3.205.055	292 / Filter discs.	100
70 mm	FTR 3.205.070	292 / Filter discs.	100
90 mm	FTR 3.205.090	292 / Filter discs.	100
110 mm	FTR 3.205.110	292 / Filter discs.	100
125 mm	FTR 3.205.125	292 / Filter discs.	100
150 mm	FTR 3.205.150	292 / Filter discs.	100
185 mm	FTR 3.205.185	292 / Filter discs.	100
240 mm	FTR 3.205.240	292 / Filter discs.	100

Grade 292a - slowly filtering, narrow pores

filtration of fine precipitates such as lead dioxide, calcium fluoride, nickel sulphide, zinc sulphide

Basis weight 97 g/m²
Filtration velocity: 10 ml / 60 seconds
Typical pore size: 5 - 8 µm
Similar to: Whatman 2, S&S 2095

Diam.	Code	Grade/Description	Pack of
55 mm	FTR 3.215.055	292a / Filter discs.	100
70 mm	FTR 3.215.070	292a / Filter discs.	100
90 mm	FTR 3.215.090	292a / Filter discs.	100
110 mm	FTR 3.215.110	292a / Filter discs.	100
125 mm	FTR 3.215.125	292a / Filter discs.	100
150 mm	FTR 3.215.150	292a / Filter discs.	100
185 mm	FTR 3.215.185	292a / Filter discs.	100
240 mm	FTR 3.215.240	292a / Filter discs.	100

GENERAL PURPOSE, QUALITATIVE-TECHNICAL FILTERS

Grade 3 hw - medium rate of filtration, increased wet strength, plain

This paper grade is generally used for analytical and technical purposes, as clarification, analytical sample preparation, particle analysis

Basis weight 65 g/m²
Filtration velocity: 10 ml / 20 seconds
Typical pore size: 8 - 12 µm
Similar to: Whatman Student Grade, S&S 595, 3002
M-N 615


Diam.	Code	Grade/Description	Pack of
55 mm	FTR 3.303.055	3 hw / Filter discs.	100
70 mm	FTR 3.303.070	3 hw / Filter discs.	100
90 mm	FTR 3.303.090	3 hw / Filter discs.	100
110 mm	FTR 3.303.110	3 hw / Filter discs.	100
125 mm	FTR 3.303.125	3 hw / Filter discs.	100
150 mm	FTR 3.303.150	3 hw / Filter discs.	100
185 mm	FTR 3.303.185	3 hw / Filter discs.	100
240 mm	FTR 3.303.240	3 hw / Filter discs.	100

GLASS MICROFIBRE FILTERS

Grade MGA

Highly efficient for general laboratory filtration, clarification of buffer and reagent solutions corresponds to many international standards for air and water pollution monitoring

Basis weight	52 g/m ²
Thickness	0,25 mm
Particle retention in liquids:	1,6 µm
Pressure drop:	33 mbar (a=10 m2, flow velocity 40 cm/sec.)
Finish:	binder-free
Similar to:	Whatman GF/A, S&S GF 50, M-N GF 1



MICRO-GLASS FIBRE FILTERS BINDER FREE UNDER A RASTER ELECTRON MICROSCOPE

Diam.	Code	Grade/Description	Pack of
25 mm	FTR 3.1101.025	MGA / Glass Microfibre discs.	100
47 mm	FTR 3.1101.047	MGA / Glass Microfibre discs.	100
50 mm	FTR 3.1101.050	MGA / Glass Microfibre discs.	100
55 mm	FTR 3.1101.055	MGA / Glass Microfibre discs.	100
70 mm	FTR 3.1101.070	MGA / Glass Microfibre discs.	100
90 mm	FTR 3.1101.090	MGA / Glass Microfibre discs.	100
100 mm	FTR 3.1101.100	MGA / Glass Microfibre discs.	100
110 mm	FTR 3.1101.110	MGA / Glass Microfibre discs.	100
125 mm	FTR 3.1101.125	MGA / Glass Microfibre discs.	100
150 mm	FTR 3.1101.150	MGA / Glass Microfibre discs.	100
180 mm	FTR 3.1101.180	MGA / Glass Microfibre discs.	100

Grade MGC

Standard filter for the clarification and monitoring of waste water and water, scintillation counting on the filter, cell harvesting, hydrocarbon analysis where cellulose fibres a nuisance as by DIN EN 872

Basis weight	52 g/m ²
Thickness	0,26 mm
Particle retention in liquids:	1,2 µm
Pressure drop:	54 mbar (a=10 m2, flow velocity 40 cm/sec.)
Finish:	binder-free
Similar to:	Whatman GF/C, S&S GF 52, M-N GF 3



Diam.	Code	Grade/Description	Pack of
25 mm	FTR 3.1103.025	MGC / Glass Microfibre discs.	100
47 mm	FTR 3.1103.047	MGC / Glass Microfibre discs.	100
50 mm	FTR 3.1103.050	MGC / Glass Microfibre discs.	100
55 mm	FTR 3.1103.055	MGC / Glass Microfibre discs.	100
70 mm	FTR 3.1103.070	MGC / Glass Microfibre discs.	100
90 mm	FTR 3.1103.090	MGC / Glass Microfibre discs.	100
100 mm	FTR 3.1103.100	MGC / Glass Microfibre discs.	100
110 mm	FTR 3.1103.110	MGC / Glass Microfibre discs.	100
125 mm	FTR 3.1103.125	MGC / Glass Microfibre discs.	100
150 mm	FTR 3.1103.150	MGC / Glass Microfibre discs.	100
180 mm	FTR 3.1103.180	MGC / Glass Microfibre discs.	100

CELLULOSE EXTRACTION THIMBLES

for Soxhlet extraction

Grade 30

extraction thimbles made of cellulose, are perfect sample reservoirs for Soxhlet extraction in quality control, research and/or analysis. Solvent extraction from solids and semi-solids can be carried out easily and reliably. The thimbles fit Soxhlet extractors and automatic extractors.

Material	Cellulose
Shape	hollow cylinders with round bottom
Max. temp. res.:	120 °C

ID x H	Code	Grade/Description	Pack of
22 x 80 mm ¹ .	FTR 1201.022080	30/extraction thimbles	25
25 x 80 mm	FTR 1201.025080	30/extraction thimbles	25
25 x 100 mm ¹ .	FTR 1201.025100	30/extraction thimbles	25
26 x 60 mm ³ .	FTR 1201.026060	30/extraction thimbles	25
28 x 22 mm ⁴ .	FTR 1201.028022	30/extraction thimbles	25
30 x 80 mm	FTR 1201.030080	30/extraction thimbles	25
33 x 80 mm ^{2,3,5} .	FTR 1201.033080	30/extraction thimbles	25
33 x 94 mm ¹ .	FTR 1201.033094	30/extraction thimbles	25
33 x 118 mm	FTR 1201.033118	30/extraction thimbles	25
43 x 123 mm ¹ .	FTR 1201.043123	30/extraction thimbles	25

1. fits Büchi B-811
2. fits Gerhardt Soxterm Automatic
3. fits Foss Tecator System
4. fits Foss Soxtec 2050
5. fits Velp Extractors SER 148



GRIDDED MEMBRANE FILTERS

Sartorius® brand


Sterile, individually packed membrane filters have long become standard for routine microbiological quality control because of the user benefits they offer. They are ready-to-use and save preparatory time, they avoid the possibility of contamination of remaining filters in opened packs and are GLP conform, having filter identification and lot number printed on each individual envelope.

Material:	Cellulose nitrate, white with black grid
Reaction to water:	Hydrophilic
Grid size:	3,1 x 3,1 mm
Pore size:	0,45 micron
Diameter:	47 mm
Presentation:	Sterile, individually packed

Diam.	Code	Pack of
47 mm	SAR 11406-47ACR	1000 (10 x 100)

Material:	Cellulose nitrate, white with black grid
Reaction to water:	Hydrophilic
Grid size:	3,1 x 3.1 mm
Pore size:	0,2 micron
Diameter:	47 mm
Presentation:	Sterile, individually packed

Diam.	Code	Pack of
47 mm	SAR 11407-47ACR	1000 (10 x 100)



SYRINGE FILTERS

Sartorius® brand

Cellulose acetate syringe filters are used for the filtration of volumes up to 100 ml. The syringe filters units are ideal for applications requiring maximum sample recovery such as: Tissue culture preparation, sterile filtration and clarification of biological fluids, specimen solutions, protein and enzyme filtration, hybridization buffers and other aqueous solutions.

Specification:

Pore size	0,20 µm	0,45 µm
Color code	blue	yellow
Membrane	Cellulose Acetate	
Housing	Cylolite	
Diameter	28 mm	
Filter aerea	6,2 cm ²	
Idle Volume	0,1 ml	
Bubble Point	3,2 bar	1,6 bar
Max. Pressure	4,5 bar	
Burst Strength	6 bar	
Max. temperature	50 °C	
Flow Rate (water)	60 ml/min	180 ml/min
Hydrophilic	yes	
Inlet	Luer lock	
Outlet	Luer lock	
Presentation	sterile, individually packed	

Pore size	Code	Color Code	Pack of
0,20 µm	SAR 16534-K	blue	50
0,45 µm	SAR 16555-K	yellow	50



ALPHABETICAL INDEX

ABBE Refractometers.....	106	Erlenmeyer Flasks with conical ground socket.....	127
Alcoholometers.....	146	Erlenmeyer Flasks, glass.....	110
Allihn Condenser.....	125	Erlenmeyer Flasks, plastic.....	136
Analytical Balances.....	18-20	Evaporation Dish, porcelain.....	131
Aspirator with Trap Flask.....	59	Extraction Thimbles.....	159
ASTM Hydrometers.....	145		
Autoclave.....	49	Filter Crucibles.....	130
Automatic Burettes, Dr. Schilling.....	123	Filter Paper.....	153-160
Automatic Burettes, Pellet.....	124	Filtering Flasks.....	110
Automatic Pipettes.....	4-7	Flat Bottom Flasks.....	110
		Flat Bottom Flasks, with conical ground socket.....	127
Balances.....	18-23	Funnels Buechner, porcelain.....	130
Baumé Hydrometers.....	145	Funnels, plastic.....	137
Beaker, glass.....	109	Funnels, short stem, glass.....	111
Beaker, plastic.....	136		
Biological Microscopes.....	96-103	Gas Burner.....	50
Blood Counting Chambers.....	114	Glass Microfibre Filters.....	158
Blood Diluting Pipettes.....	114	Gooch Crucibles, porcelain.....	129
Bottle-Top Burette, digital.....	16	GP Series Micropipettes.....	4-5
Bottle-Top Dispenser.....	14-15	Graduated Cylinders, glass.....	118
Bottles, weighing.....	126-127	Graduated Cylinders, plastic.....	134
Brix Hydrometers.....	145	Graduated Pipettes.....	114-115
Buechner Funnels, porcelain.....	130	Graduated Pitchers, plastic.....	137
Buffer Solutions.....	76		
Burettes, Daffert.....	123	Haematology Pipettes.....	114
Burettes, Mohr.....	116, 122	Hand Refractometers.....	105
		Heating Mantles.....	58
Centrifuge Tube, conical, with screw cap.....	142	Hellendahl, Staining Jar.....	113
Centrifuges.....	80-94	Haematocrit Centrifuge.....	80
Chemical Resistant Vacuum Pumps.....	61	Hydrometers.....	144-147
Clinical Photometer.....	31		
CO2 Incubator.....	48	Incubator.....	43
Coiled Condenser.....	125	Inverted Biological Microscope.....	97
Colony Counter.....	63		
Colour Comparison Tubes.....	119	Jacket Coil Condenser.....	126
Combispin, Centrifuge & Vortex.....	37		
Condenser, glass.....	125-126	Laboratory Bottles, ISO, clear & amber.....	109
Conductivity/TDS/Salinity/Temp. Meters.....	78-79	LAC Milk Analyzer.....	34-35
Conductivity Calibration Solutions.....	76	Lactometers.....	146
Conical Test Tubes.....	141	Lids for porcelain crucibles.....	129
Coplin, Staining Jar.....	113	Liebig Condenser.....	125
Counting Chambers.....	114		
Crucibles, porcelain.....	129-130	Magnetic Filter Holders.....	62
Cuvettes.....	30	Magnetic Stirrer.....	39
Cylinders, graduated, glass.....	118	Magnetic Stirrer with Hotplate.....	39-40
Cylinders, mixing.....	119	Magnetic Stirring Bars & Retrievers.....	139
		Manifolds, stainless steel.....	62
Density Hydrometers.....	144	Max/Min Thermometers.....	150
Desiccator Plates, porcelain.....	111, 130	Measuring Cylinders, glass.....	118
Desiccators.....	111	Measuring Cylinders, plastic.....	134
Digital ABBE Refractometer.....	106	Measuring Pipettes.....	115
Digital Bottle-Top Burettes.....	16	Membrane filter.....	159
Digital Hand Refractometers.....	105	Microburettes, Bang.....	123
Digital Microscope Cameras.....	104	Microliter Centrifuges.....	86-88
Dimroth Condenser.....	126	Microfibre Filters, glass.....	158
Dispensers, Bottle-Top.....	14-15	Micropipettes.....	4-7
Disposable Cuvettes.....	28, 142	Microplate Incubator Shaker.....	43
Disposable Serological Pipettes, plastic.....	142	Microplate Shaker.....	43
DO/Temp. Meter.....	77	Microplate reader.....	32
Dropping Bottles, glass.....	112	Microplate washer.....	33
Dropping Funnels.....	121	Microscope Cameras, digital.....	104
Dry Bath Block Incubator.....	54	Microscopes.....	96-103
		Microtest Plates.....	140
Electrolyte Solution.....	73	Microtest Tubes.....	140
Electronic Pipette Controller.....	13	Milk Analyzer.....	34-35
Electronic Pipettes.....	8	Mini-Rocker Shaker.....	42

ALPHABETICAL INDEX

Mini-Shaker, Sunflower type	38	Series Pipettor	9
Mixing Cylinders, glass	119	Shaker for Microtubes, with heating & cooling	55
Moisture Analyzer	22	Soxhlet Extractor	126
Monocular Microscopes	102	Soxhlet Extractor Apparatus	126
Mortars, porcelain	132	Spectrophotometers	25-27
Multi Rotator	42	Square Bottles, plastic	138
Multistand Rack	138	Staining Dishes	113
Muffle Furnace	52	Stereo Microscopes	103
		Stereo Zoom Microscopes	103
Nessler Tubes	119	Stoppers, glass	127
Narrow Neck Bottles, plastic	138	Stoppers, plastic	141
		Stopwatch	148-149
Optical and Quartz Cuvettes	30	Sunflower Mini-Shaker	38
Orbital Shaker	44	Syringe Filter	160
Orbital Shaker Incubator	46		
Overhead Stirrer	41	Test Tube Racks	138
		Test Tubes	114-115
Pasteur Pipettes, glass	114	Test Tubes with thread	115
Pasteur Pipettes, plastic	141	Test Tubes, plastic	140
PCR Tubes	140	Thermal Cycler	56-57
PCR Workstation	53	Thermo Shaker for Microtubes	55
PCR Thermo Cycler	56-57	Thermo-Alcoholometers	147
Pestles for Mortars, porcelain	132	Thermo-Hygrometers, digital	152
Petri Dishes, glass	112	Thermo-Lactometers	146
Petri Dishes, plastic	140	Thermometers, digital	150
pH Electrodes	73	Timer, mechanical & digital	148-149
pH Indicators	77	Tube roller	38
pH Meters	73		
Photometer	28	Universal Centrifuges	80-94
Pipette Aids	12	Universal Orbital Shaker	44
Pipette Controller	12	UV/VIS Spectrophotometer	24-27
Pipette Controller, electric	13		
Pipette Fillers	12	Vacuum Pumps	61
Pipette Rack	139	Vials with thread	115
Pipette Tips	10-11	VIS Spectrophotometer	24-25
Pipettes, glass	114, 117	Volumetric Flasks for sugar testing	121
Pipettes, plastic	141	Volumetric Flasks, glass	120
Pipettes, Haematology	114	Volumetric Flasks, plastic	135
Pitchers, graduated, plastic	137	Volumetric Pipettes	115, 117
Pouring Rings	109	Vortex Mixer	36
Precision Balances	18-21		
Programmable Orbital Shaker	45	Wash Bottles	138
Pycnometers	117	Water Bath	51
		Water Deionizer	72
Qualitative Filter Paper	156-157	Water Purification System	64-69
Quantitative Filter Paper	154-155	Water Stills	70-72
		Watch glasses, glass	112
Rack for Pipettes	139	Weighing Bottles	126-127
Rack for Test Tubes	138	Westergreen Pipettes	114
Reagent Bottles, glass, clear & amber	112	Wide neck bottles, plastic	138
Red Rubber Bulbs for Automatic Burettes	125	Wintrobe Pipettes	114
Refractometers	105-106		
Reservoir Bottles, for Automatic Burettes	125		
Rocker Shaker	42		
Rotary Evaporator	60		
Rotator	37, 42		
Round Bottom Flasks	110		
Round Bottom Flasks, with conical ground socket	128		
SA Series Micropipettes	6-7		
Safety Laboratory Gas Burner	50		
Sahli Pipettes	114		
Schieferdecker, Staining Dish	113		
Screw caps	109, 115		
Sedimentation Cones	119		
Sedimentation Cone, plastic	137		
Separating Funnels	121		